IFU-OTH-0022

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Revisão / Revision: 1
Data / Date: 30/01/2026

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Technical Responsible: Dr. Silvio Zilio - CRO-SP n°: 47144

Product description

The Nano Tube MJanson is a bonding accessory, manufactured in polycrystalline ceramic (Al₂O₃) using CIM (Ceramic Injection Molding) technology and features a base with a design that consists of a base with pins without surface treatment, with double curvature, being mesiodistal and cervical-occlusal, providing good adaptation to the tooth surface and the presence of pins guarantees good adhesion. It was designed for use with arches up to 0.018". This accessory is 3.50 mm wide and 2.20 mm high. It is characterized by rounded corners and translucency, resulting in comfort and aesthetics for the patient. Manufacturing technology offers professionals quality, safety and reliability when carrying out clinical procedures.

The Nano Tubo is an accessory designed to be used as an aid in mechanics with aligners and small orthodontic movements with fixed appliances.

**Illustrative images*

Intended use

The Nano Tube is used in hybrid mechanics with aligners in positions previously determined by the professional according to the mechanics to be performed and, according to treatment planning, acting to establish torque with the use of round nickel-titanium arches with a section of up to 0.018". They can also be used to correct rotations, crowding recurrences and alignment and leveling of molars with short or semi-erupted clinical crowns. They can be used in patients with:

I- Sagittal problems:

- 1) Angle Class I malocclusion;
- 2) Angle Class II malocclusion;

3) Angle Class III malocclusion;

II- Inter-arch and intra-arch problems:

1) Dental crowding;

Indication for use

The Nano Tubo is used in one or more dental units, together with aligners or alone, in positions determined by the professional according to the mechanics to be used and according to treatment planning. They are generally used together with round nickel-titanium archwires with a section of up to 0.018", according to the intended use mentioned above.

Clinical indications - Deviations from ideal occlusion (mentioned in Intended Use) and the associated psychosocial issues today play an important role in defining the need for orthodontic treatment.

How to use:

For proper bonding of the Nano Tube, the following sequence of procedures must be adopted and are summarized as:

1. prophylaxis;
2. acid etching of the enamel;
3. application of the primer followed by bonding the Nano Tube to the tooth.

The instructions for bonding are described below:

I-Direct bonding

- 1) Isolation of the work field with plastic cheek retractors and relative isolation with cotton rollers, maintaining the high-power saliva sucker, to avoid contamination of the bonding areas with saliva;
- 2) Delimitation of the specific area for gluing the accessory;
- 3) Carry out prophylaxis with pumice stone, water and Robson's brush to remove the Nasmyth membrane from the enamel surface; the acquired film and other contaminating agents;
- 4) Wash the sanitized areas with plenty of water;
- 5) Condition the enamel surface, only the area that will receive the Nano Tube, using 37% phosphoric acid gel applied for 15 to 30 seconds;
- 6) Wash the conditioned areas intensively, keeping the high-power sucker in position, to avoid contamination of these enamel surfaces. If contamination with saliva occurs, the conditioning procedure is repeated;
- 7) Drying the conditioned enamel surface with a jet of air free from contamination with water or oil;
- 8) Application of fluid resin or primers to the enamel surface in relation to the bonding area with a brush or microbrush, without exerting pressure;
- 9) Photopolymerization with the light-curing device, at a time and ultraviolet emission defined by the professional, according to the treatment;

- 10) Keep the accessory fixed between the ends of the tweezer for Nano Tube, avoiding manipulating the accessory directly with your fingers, which leads to contamination;
- 11) Application of a thin layer of primers over the entire base of the Nano Tube, covering the entire retention area, being an optional procedure;
- 12) Positioning the Nano Tube on the enamel surface using the tweezer for Nano Tube, applying light pressure;
- 13) Removal of excess resin; and
- 14) Light-curing with the LEDX-T 2400 light-curing device, in the fast Ortho function, for 3 seconds per button or another light-curing device with similar function and power. An alternative light-curing device that has similar technical characteristics and is properly calibrated can be used.

Due to the large number of steps involved in the accessory bonding process, it becomes very difficult to control the moisture absorbed in the process, which may lead to contamination of the enamel surface and, consequently, failure of adhesion to the tooth. The presence of moisture and/or contamination during the bonding process is the biggest cause of adhesion failures, delaying treatment and increasing costs.

The Nano Tube may undergo changes in its surface structure when exposed to adverse conditions present in the oral environment, especially foods with dye and pigments. Therefore, the patient must be guided regarding correct hygiene and be aware of the risks of staining if adequate care is not taken.

The medical product in question can only be placed on a patient by a dental surgeon specialist in Orthodontics.

Nano Tube Removal

After achieving the proposed treatment objectives, the accessory is removed, however, without damaging the tooth enamel structure. It is recommended to use appropriate and good quality instruments to avoid discomfort to the patient and the risk of injury to soft tissues and dental structures.

The technical removal procedures involve the following steps:

- 1) Removal of excess resin that may be surrounding the base of the Nano Tube, however, if the bonding has been carried out with technical care, the minimum possible amount of resinous material around the orthodontic accessory is observed;
- 2) Use of How Reto (number 110) or Weingart (number 120) pliers, both used to perform mesiodistal compression of the base of the accessory, which promotes deformation of the accessory and release of it from the enamel surface, through the rupture in the Nano Tube-resin interface. The instruments used at this stage will depend on the professional's familiarity, with some using 346 pliers or metallic bracket puller pliers. Removal must be carried out with control and precision to avoid the risk of swallowing or aspiration.
- 3) Remove the adhesive using scraper pliers with Widia or titanium tips. Excess resin can also be removed with a pair of pliers;
- 4) Removal of composite residues with a multilaminar conical trunk drill made of tungsten carbide, adapted to a 9-blade low rotation contra-angle. It can also be used the high-rotation 24-blade drill, code CF 375R or similar of the same quality. These drills have inverted blades to preserve tooth enamel; and
- 5) After removing the resin residue, the surface of the enamel should be finished with a Shofu Ceramisté polishing rubber. The final polishing is carried out with rubber cups and Spain white in water or cloth wheels in a contra-angle and diamond paste.

Composition:

Nano Tube – Aluminum oxide

Storage Conditions:

In intact and closed packaging.

This product does not require temperature, humidity, or light control during storage. The recommendation is to keep it in a clean, dry place.

The manufacturer recommends the single use of this product. It should be disposed of correctly in accordance with the information contained in the Manual of Dental Services Risk Prevention and Control, Technology in Health Services, 2006 (available on the ANVISA website).

Warnings and Symbology:

- It is essential that only the specialist installs the product on the patient, since in addition to the clinical examination, the professional must have the necessary knowledge to make a systematic description of the patient's dentofacial features, list the priority problems, describe possible solutions for each of the problems, considering the possible favorable or unfavorable interactions for each of the proposed attempts.
- It is also the orthodontist's duty to determine the maximum expected gain from orthodontic therapy and the minimum loss, or to calculate the benefit versus risk ratio for each patient, with both risk and benefit including psychosomatic factors to determine the best management for each patient.
- Patients should be instructed in oral hygiene, maintenance, and food care. The following should be avoided: ice, coconut, hard and/or sticky sweets, caramels, nougats, jellybeans, corn, popcorn and sugar cane. Biting hard objects such as pens and pencils has the same effect. Chewing gum and toffee often loosen the bands or rings that have been cemented in place.
- Advise patients NOT to chew or bite on hard foods during orthodontic treatment to prevent premature detachment. This product should not be used in an MRI environment, as it may cause interference in the images during the examination. It is up to the orthodontic professional to warn the patient of this contraindication and arrange for the appliance to be removed if necessary.
- The high concentration of bacteria and fungi present in the mouth can accelerate the corrosion of metal orthodontic appliances. The pH of the oral environment also has a major influence on the corrosive process, so a diet rich in acidic foods, such as excessive intake of soft drinks with acidic pHs, should be avoided as much as possible during treatment.
- Do not touch the Nano Tube with ultrasonic tips: this can lead to fractures and detachment.
- The product is supplied in non-sterile condition, however it has proven biocompatibility and no cleaning, disinfection or sterilization measures are required, provided the packaging has not been tampered with. If the product comes into contact with surfaces or substances other than those intended for use, it is recommended that it be discarded.
- The product must not be reused or reprocessed. Do not use the product if its packaging is damaged or open, to minimize possible dangers related to the transmission of microorganisms or possible product counterfeiting.

- It is the responsibility of qualified professionals to advise their patients of the risk of aspiration and swallowing in the event of breakage/ detachment of parts, not to underestimate the possible immediate or delayed complications resulting from this type of accident and to adequately inform the patient and/or their legal guardians of the risks inherent in the event and the measures to be taken.
- This product should only be used by qualified orthodontic professionals.
- Use by people without the necessary technical knowledge could lead to unwanted tooth movement, bone loss, tooth loss, fenestration, gingival recession, or root dehiscence.
- It is the responsibility of the qualified professional to advise the patient on correct oral hygiene, in order to avoid the appearance of plaque and tartar or diseases such as gingival inflammation (gingivitis), periodontitis or even endocarditis. Considering that the product is used in contact with mucous membranes and body fluids, it is recommended that professionals apply good clinical practices for the disposal of medical products in accordance with current health regulations.

ATTENTION: Orthometric cannot be held responsible for the results obtained by applying the product without the necessary precautions or failing to observe the warnings.

Contraindications

Orthometric only supplies its products to qualified professionals. It is the orthodontist's sole responsibility to identify conditions that may be contraindicated for treatment, such as:

- Patients with poor oral hygiene.
- Patients unable to cooperate with treatment.
- Patients with an oral environment conducive to the demineralization of tooth enamel.

	Do not re-use		Consult instructions for use
	Keep away from sunlight		Do not use if package is damaged
	Medical device		Caution
	Catalogue number		Non-sterile

IFU-OTH-0022

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Revisão / Revision: 1
Data / Date: 30/01/2026

	Batch code		Keep dry
	Manufacturer		Date of manufacture
	Use-by date		