Pt

Es

En

Fr

Ethestic Plus

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Ethestic Plus

Silicone brushes for shaping composites

Ethestic Plus were developed to promote easiness and precision in shaping composites, both direct and indirect restorations. The kit contains a straight instrument, which facilitates the access to the vestibular anterior/posterior teeth, as well as for laboratory work, and an angled instrument, which allows an easier and more comfortable access to the lingual/palatal region of the anterior/posterior teeth and to the occlusal region of posterior teeth. The brushes are interchangeable, made of silicone rubber, anti-heat and available in two different shapes and two different hardness (hard and soft tops), which allows to shape the composite to obtain the desired tooth anatomy. The possibility of rotating the brushes on the end of the instrument provides comfort during use. Using SepticCom Shape avoids the typical lines and marks caused by hard instruments and, consequently, the differences of shine and color that they may cause.

Use: Select the instrument and brush which best fits to the clinical case. Place the brush on the end of the instrument. Rotate the brush to reach the best position for the work. The composite with light touches, starting with the hard brushes (blue) and finishing with the soft brushes (white).

Warnings, precautions and restrictions: Professional use. If there is proven allergy to any of its components, the product must not be used. To avoid an accidental swallowing of this product, it is necessary the application of a total isolation technique of the operatory and it is recommended to clean the silicone brushes with alcohol after each use to remove composites excess in order to ensure their functionality. Product not recommended for condensable composites, flow composites, amalgams, glass ionomer cement and self-curing acrylic resins. Any serious incident related to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Cleaning and sterilization: Autoclavable (126°C - 16min / 134°C - 18min). Before sterilization, proceed with the cleaning using water and enzymatic detergent. Cleaning: Carry out the cleaning with the aid of water and enzymatic detergent 4 enzymes (Amylase, Protease, Lipase and Carboxydase) and with pure pH range between 6.5 - 7.5. Prepare the solution according to the manufacturer's instructions. The water should be potable. For a better effectiveness, it is recommended that the temperature be between 40°C and 65°C. Submerge the products and leave in contact for 5 minutes with the solution. Remove parts and rinse under running water for at least 10 seconds. Dry the pieces with absorbent paper or disposable wipes. Inspect: Check the product before each use, seeing if there is damage to brushes that prevent proper performance, and do not use in such cases.

Storage: To be stored in a cool and single place, protected from sunlight.

Dimensions: Available in 4 different shapes: double chamfered edge, double angled chamfered edge, cut, double angled chamfered edge, cut, flange.

Composition: Instruments: Polysulfone. Soft brushes: Silicone and titanium dioxide. Hard brushes: Silicone, titanium dioxide, aluminum oxide and atoxic colorants.

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