200 mm

Liner

PORTUGUÊS

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ENGLISH

1. Product Description
Liner is a self-cure, rigid and radiopaque compound based on calcium hydroxide, indicated for pulp and protective lining under restorative materials, cements and other materials.
Liner does not inhibit the polymerization of acrylic and composite restoration. It has high resistance to dissolution by phosphoric acid, acid conditioning technique and optimum solubility in water.

1.1. Indications
Liner is indicated for pulp capping and protective lining under restorative materials, cements and other materials.

1.2. Composition
Base Paste: Barium Sulfate, Silicon Dioxide, Titanium Dioxide, 1,3-Butylene Glycol Disalicylate and Pigments.
Catalyst Paste: N-ethyl-p/o-toluene sulfonamide, Calcium Hydroxide, Zinc Oxide, Zinc Stearate and Pigment.

1.3. Presentation
1tube of base with 13g
1tube of catalyst with 11g
1Instructions for use
1mixing pad

2. Instructions for use
1 - Dispense equal amounts of base And catalyst (1: 1) on the block of mixture.
2 - Immediately mix the base with the Catalyst, using an instrument of The active tip rounded to a Uniform color. Do not mix overly. Usually one time From 10 to 18 seconds are sufficient for Obtaining a homogeneous mixture.
Note: Working and prey times Vary with changes in temperature And humidity. With increasing humidity and Of the temperature are reduced the time of Available work and

the holding time. On the contrary, the reduction of humidity and Reduce the temperature increase The available working time and the Prey Under ambient conditions Normal, the material will take 2'30 "and 3'00" minutes on the block of mixture. Hold time is shorter in the mouth.
3- Apply the mixed material inside The cavity, dispensing it with Circular movements on walls Cavity, obtaining a thickness of Minimum of 0.75 to 1.5 mm of cement. THE Cavity should be dry.
4- After the prey, remove any Excess material from the areas of Retention of the enamel margins as Also of the cavo-superficial angle.

3. Precautions, Contraindications and Warnings

3.1. Precautions

- All dental material must only be used by dentists or authorized dental laboratory technician.
- Keep the product away from children.
- Keep tubes tightly closed.
- Avoid contact with skin and eyes. If accidental of contact with skin occurs, wash immediately with plenty of soap and water. In case of contact with eyes, rinse immediately with plenty of water and seek medical attention.
- May cause allergic reactions in persons susceptible to sulfonamides.
- Acetone is detrimental to calcium hydroxide materials. Minimize the contact with acetone.
- Minimize contact of these products with calcium hydroxide liners.

3.2. Contraindications and Warnings
There are no restrictions when Liner is used by a certified professional following the Directions for Use. Avoid using Liner on any patient with allergic history to sulfonamides or other components of the product.

4. Validity
Do not use after expiration date. If stored appropriately, Liner has a shelf life of 36 months.

5. Storage
Stored at room temperature (15°C - 30°C) or below in a dry place and protected from light and heat. Separation of solid and liquid phases of the base paste is accelerated if the material is exposed to

high temperatures. When stored at reduced temperatures, allow the material to return to normal room temperature before using.

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ENGLISH

1. Product Description

Liner is a self-cure, rigid and radiopaque compound based on calcium hydroxide, indicated for pulp and protective lining under restorative materials, cements and other materials.

Liner does not inhibit the polymerization of acrylic and composite restoration. It has high resistance to dissolution by phosphoric acid, acid conditioning technique and optimum solubility in water.

1.1. Indications

Liner is indicated for pulp capping and protective lining under restorative materials, cements and other materials.

1.2. Composition

Base Paste: Barium Sulfate, Silicon Dioxide, Titanium Dioxide, 1,3-Butylene Glycol Disalicylate and Pigments.

Catalyst Paste: N-ethyl-p/o-toluene sulphonamide, Calcium Hydroxide, Zinc Oxide, Zinc Stearate and Pigment.

1.3. Presentation

1tube of base with 13g

1tube of catalyst with 11g

1Instructions for use

1mixing pad

2. Instructions for use

1 - Dispense equal amounts of base And catalyst (1:1) on the block of mixture.

2 - Immediately mix the base with the Catalyst, using an instrument of The active tip rounded to a Uniform color. Do not mix overly. Usually one time. From 10 to 18 seconds are sufficient for Obtaining a homogeneous mixture.

Note: Working and prey times Vary with changes in temperature And humidity. With increasing humidity and Of the temperature are reduced the time of Available work and

the holding time. On the contrary, the reduction of humidity and Reduce the temperature increase The available working time and the Prey Under ambient conditions Normal, the material will take 2'30 "and 3'00" minutes on the block of mixture. Hold time is shorter in the mouth.
3- Apply the mixed material inside The cavity, dispensing it with Circular movements on walls. Cavity, obtaining a thickness of Minimum of 0.75 to 1.5 mm of cement. THE Cavity should be dry.

4- After the prey, remove any Excess material from the areas of Retention of the enamel margins as Also of the cavo-superficial angle.

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3.1. Precautions

• All dental material must only be used by dentists or authorized dental laboratory technician.

• Keep the product away from children.

• Keep tubes tightly closed.
• Avoid contact with skin and eyes. If accidental of contact with skin occurs, wash immediately with plenty of soap and water. In case of contact with eyes, rinse immediately with plenty of water and seek medical attention.

• May cause allergic reactions in persons susceptible to sulfonamides.

• Acetone is detrimental to calcium hydroxide materials. Minimize the contact with acetone.

• Minimize contact of these products with calcium hydroxide liners.

3.2. Contraindications and Warnings

There are no restrictions when Liner is used by a certified professional following the Directions for Use. Avoid using Liner on any patient with allergic history to sulfonamides or other components of the product.

4. Life

Do not use after expiration date. If stored appropriately, Liner has a shelf life of 36 months.

5. Storage

Stored at room temperature (15°C - 30°C) or below in a dry place and protected from light and heat. Separation of solid and liquid phases of the base paste is accelerated if the material is exposed to

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