It has in its composition a monomeric organic matrix whose polymerization process occurs through the formation of free radicals. The camphorquinone present in the composition is activated by the incidence of visible light with a wavelength varying between 400 and 500 nm (blue spectrum) with maximum absorption at 468 nm. When activated, this component transfers one of its electrons to a supporting substance in the photoinitiation, resulting in the formation of free radicals and triggering the polymerization process. Several factors contribute to an efficient light-curing: incidence of light in the spectrum from 400 to 500 nm, power of the light-curing unit, and distance from the light source, exposure time and thickness of the material layer.

The curing time is 40 ± 10 seconds if a light-curing unit is used with a wavelength of 400 to 500 nm, with a maximum emission peak of 468nm (Blue Light) and a minimum power of 400mW/cm^2 , totaling a supply of at least 16J of energy if the curing distance is less than 5 mm.

USE INSTRUCTIONS:

I - DEFINITIVE AESTHETIC DENTAL RESTORATIONS:

1. Select colors with the moistened teeth;
2. Promote coronary polishing with pumice and water;
3. Perform the cavity preparation according to the recommended technique;
4. Perform absolute isolation;
5. Perform the etching (e.g. ATTAQUE GEL) of the enamel for 30 seconds and dentin for 15 seconds;
6. Wash the etched area with plenty of water for 30 seconds;
7. Dry the dentin (without dehydrating) with light intermittent jets of air. It is important to dry the enamel well and not to dehydrate the dentin;
8. Use the adhesive system (e.g. MASTER BOND), following the manufacturer's guidelines;
9. Apply the resin in thin layers of 1 to 2 mm, light curing each layer;
10. Promote finishing and polishing;
11. Perform occlusion control;
12. Polish with sanding discs and abrasive silicone tips;
13. It is recommended to apply surface sealant for restorations (e.g. BIOFORTY), following the manufacturer's guidelines.

II - REPAIR OF PROSTHESES AND INDIRECT RESTORATIONS:

1. To increase adhesion, apply the porcelain conditioner, for the time recommended by the manufacturer, wash and dry;
2. Apply SILANO over the fractured porcelain and wait for 1 minute; perform a new application;
3. Use the adhesive system (e.g. MASTER BOND), following the manufacturer's guidelines;
4. Apply the resin in thin layers of 1 to 2 mm;
5. Promote finishing and polishing.

* Perform, if necessary, the masking of the metal.

SPECIAL CARES:

- The dentin-pulp complex must be protected with an appropriate base before etching.
- Eugenol-based tools have a negative influence on the light-curing time.
- The use of cationic mouthwashes, plaque detectors and chlorhexidine can interfere with the color of the material.
- For satisfactory light-curing, use a light-curing unit of adequate quality, and the distance for the light-curing of the restoration surface must be less than 5 mm; carefully observe the exposure time, as well as the material thickness; Resins not properly cured can cause irritation to the dentin-pulp complex, as well as reducing mechanical properties and causing color change. In long-term work, deflect the light from the reflector to avoid premature hardening of the material.
- It is recommended to use absolute insulation.
- It is recommended to apply surface sealant for restorations (e.g. BIOFORTY), to seal the margins of restorations and increase the longevity of restorations.
- Dispose of the product in accordance with local legislation, ensuring its complete lack of characterization in order to prevent its reuse and environmental damage.
- Dispose the product in accordance with local regulations, ensuring its complete distortion, preventing its reuse and environment damages.

CONTRAINDICATION:

EPIC is contraindicated for people with reported sensitivity to any of the product's components.


WARNINGS:

Do not use material after expiry date.

The manufacturer is not liable for damages caused by improper use or use not provided at these instructions.

"KEEP MATERIAL OUT OF THE REACH OF CHILDREN".

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USO EXCLUSIVO DEL CIRURGÃO DENTISTA

EPIC

Light-curing Nano-hybrid Composite Resin

biodinâmica

CONTENTS:

1X 4g syringe Epic Body / Universal in colors: A1; A2; A3; A3,5; A4; B1; B2; B3; B4; C1 and C2.

1X 4g syringe Epic Enamel in colors: EA1; EA2; EA3; EA3,5; EA4; EB1; EB2; EB3; EB4; EC1; EC2; EC3; EC4; ED2; ED3 and ED4.

1X 4g syringe Epic Dentin in colors: DA1; DA2; DA3; DB1; DB2 and Dc2.

1X 4g syringe Epic Effects: Translucent; Incisal; WB and WO.

1X 2g syringe pigment in colors: Yellow, Blue, White and Brown.

Reg. ANVISA: 10298550210

COMPOSITION:

Matrix of methacrylic monomers, inorganic fillers, initiators, stabilizer, preservative and pigments.

INDICATION:

EPIC is a nano-hybrid resin indicated for definitive esthetic dental restorations anterior and posterior; fixing teeth with mobility, prosthesis repairs and indirect restorations.

TECHNICAL INFORMATION:

EPIC is a high-density, light-curing, radiopaque and fluorescence nano-hybrid composite resin. It contains in its inorganic composition a mixture of micro and nanoparticles, ranging in size from 5 nm to 2 μm. This size distribution allows for more efficient charge compaction, where the smaller particles fill the spaces between the larger particles. The inorganic composition is equivalent to approximately 79% by weight.

EPIC has fluorescence and excellent finishing and polishing, as well as excellent mechanical properties, which can be used both for previous restorations, including facets, and for posterior restorations. Complete color system, with colors for enamel, dentin, body, effects and pigments.