4. Introduce the endodontic cement into the root canal according to the selected endodontic technique.  
Note: If a more fluid cement consistency is required to facilitate insertion into the root canals, heat the glass plate at 10 to 15 cm from a flame. This procedure can be repeated as many times as necessary.

#### SPECIAL CARES:

- This product may change colour over time. It is important to clean the pulp chamber with alcohol after obturation, to remove cement residue, remaining product residues in the pulp chamber may change the colour of the tooth.
- If an irrigating solution is used in the root canals, and contains peroxides, eg. hydrogen peroxide, re-irrigate with sodium hypochlorite and saline should be performed, dry the duct well before applying HYDROSEALER.
- Immediately after use, thoroughly clean instruments used in handling material with alcohol, acetone or chloroform to avoid solidification, which may make it difficult to remove material residue.
- Avoid contact with the skin. In case of contact, immediately wash with plenty of soap and water.
- In case of contact with eyes, rinse immediately with plenty of water and seek medical advice.
- Dispose of the product in accordance with local legislation, ensuring its complete decharacterization in order to prevent its reuse and environmental damage.

#### CONTRAINDICATION:

**HYDROSEALER** is contraindicated for people with reported sensitivity to any of the product's components.


#### WARNINGS:

Do not use products beyond the stated expiration date.

The manufacturer is not liable for damage caused by incorrect use or not foreseen in the instructions for use.

" KEEPING PRODUCT AWAY FROM CHILDREN REACH "

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ANVISA Reg.: 10298550130.

### COMPOSITION:

**HYDROSEALER Powder:** Bismuth Trioxide, Bismuth Subcarbonate, Calcium Hydroxide, Hexamethylene Tetramine, Titanium Dioxide.

**HYDROSEALER Liquid:** Bisphenol Epoxy.

### INDICATION:

HYDROSEALER is an endodontic cement for root canal filling.

### TECHNICAL INFORMATION:

HYDROSEALER is an endodontic cement containing calcium hydroxide as one of the constituents of the formulation, which stimulates restorative dentin deposition, promotes an bactericidal and bacteriostatic effect and keep the alkaline pH (pH +.11).

HYDROSEALER exhibits no contraction during the setting reaction, has low dimensional stability, low solubility, high adhesiveness, sealing ability, high flow, excellent physical and biological properties, with ease of work and high radiopacity.

### INSTRUCTIONS FOR USE:

- After completion of the instrumentation, clean and dry the root ducts;
- Place HYDROSEALER on a 2: 1 resin / powder ratio glass plate (1 dust meter for 2 resin drops). The dispenser that comes with the product should be used to dose the amount of powder, the flask of resin should be positioned vertically so that the liquid is expelled in the optimal amount;
- Mix the product by incorporating the powder into the resin until a smooth and homogeneous mixture is obtained. The optimal consistency is achieved when a strand is formed, approximately 1.5 to 2.5 cm long, moving the spatula away from the plate mixture;