material and the application is simple which can be done over moist dentin due to its hydrophilic property. It promotes a high strength adhesion and protect against micro-leakage.

It's essential the use of etching agent, 37% Orthophosphoric Acid (a.e. Ataque Gel).

USE INSTRUCTIONS:

The use of rubber dam is highly recommended during all procedures.

I. In Restorations with Composites:

1. Apply the 37% Orthophosphoric Acid (a.e. Ataque Gel) - 15 seconds for dentin and 30 seconds for enamel;
 2. Rinse thoroughly with water during 30 seconds;
 3. Dry dentin (no dehydrate) with absorbent paper. It is very important a complete dry of the enamel. Dentin must shine some moist glaze;
 4. Apply **EXTRA BOND DE** with a brush (uniformly) on the etched surface, apply air jet to evaporate the solvent;
 5. Then, apply another layer of **EXTRA BOND DE** such as step 4;
 6. Light cure for 20 seconds;
 7. Apply restorative material (a.e. Master Fill) according to the manufacturer's instructions;
- II. In Ceramics, Metal and Composite (Inlays and Complete Crown) Indirect Restorations:
1. Apply SILANE on the internal surface of the ceramic parts and stand still for one minute; dry with airflow; apply new layer;

2. Prophylaxis with pumice powder and water the prepared tooth;
3. Apply the 37% Orthophosphoric Acid (a.e. Ataque Gel) for 15 seconds;
4. Rinse thoroughly with water during 30 seconds;
5. Dry dentin (no dehydrate) with absorbent paper; Dentine must shine some moist glaze;
6. Apply **EXTRA BOND DE** with a brush (uniformly) on the etched surface, apply air jet to evaporate the solvent;
7. Then, apply another layer of **EXTRA BOND DE** such as step 6;
8. Light cure for 20 seconds;
9. Use composite cement according to the manufacturer's instructions.

III. In the Prepare of Metal-Ceramic Restorations:

1. Roughen the exposed metal and ceramic surface with aluminum oxide jet or diamond burs;
2. Bevel the fractured porcelain margins;
3. Rinse and dry;
4. Apply the metal opacifier (a.e. Biocal);
5. Apply the Porcelain Conditioner for 4 minutes, rinse and dry;
6. Apply Silane all over the conditioned area and wait for one minute; apply new layer;
7. Apply **EXTRA BOND DE** uniformly, all over the conditioned area, apply air jet to evaporate the solvent;
8. Apply new layer of **EXTRA BOND DE** such as step 7;
9. Light cure during 20 seconds;
10. Complete the repair with composite filling material (a.e. Master Fill) according to the manufacturer's instructions.

Note: Masking the metal if necessary.

IV. Porcelain Facets Cementing:

1. Veneer's internal surface must come from Prosthesis lab etched with Hydrofluoric acid;
2. Place Silano and wait during one minute; apply new layer;
3. Prophylaxis the dental element with pumice powder and water;
4. Rinse and dry;
5. Etch with 37% phosphoric acid (a.e. Ataque Gel) for 15 seconds;
6. Rinse thoroughly with water during 30 seconds;
7. Dry dentin (no dehydrate) with absorbent paper. It is very important a complete dry of the enamel. Dentin must shine some moist glaze;
8. Apply **EXTRA BOND DE**, immediately and uniformly, all over the conditioned area, apply air jet to evaporate the solvent;
9. Apply new layer of **EXTRA BOND DE** as indicated in step 8;

10. Light cure during 20 seconds;
11. Mix the cement and place it in the veneers;
12. Place the veneer carefully on the teeth and remove the excess;
13. Light cure firstly in the cervical margin and then the crown and incisal margins. Avoid any contact with any part of the light curer.

SPECIAL CARES:

Manipulate all material wearing latex gloves to prevent any skin contact;

If the material gets in contact with the eyes, rinse them immediately with water;

In case of discomfort, seek medical advice;

The risk of post-operative sensitivity can be minimized with the conservative cavity preparation;

The restorative material must be placed and light cured in layers to decrease the light polymerization contraction and occlusion must be always checked. During the use of light curer, follow the manufacturer's instructions regarding eyes protection;

After use keep flask tightly close to avoid exposure to light;

Do not store **EXTRA BOND DE** near products containing Eugenol. Eugenol interferes in the polymerization.

Dispose the product in accordance with local regulations, ensuring its complete distortion, preventing its reuse and environment damages.

CONTRAINDICATION:

EXTRA BOND DE should not be used in patient is known to have sensibility to any of the ingredients listed.


WARNINGS:

Do not use material after expiry date.

The manufacturer is not liable for damages caused by improper use or use not provided at these instructions.

“KEEP MATERIAL OUT OF THE REACH OF CHILDREN”.

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EXTRA BOND DE

Light Curing Single Component
Adhesive for Enamel and Dentin

biodinâmica[®]

CONTENTS:

1x 4mL, 1x5mL or 1x6mL flask **EXTRA BOND DE**

Reg. ANVISA: 10298550196.

COMPOSITION:

Bisphenol A glicidylmethacrylate (43,79%); Methacrylate Groups (14,42%); Ethil Alcohol, Sodium Fluoride and Catalysts.

INDICATION:

EXTRA BOND DE is used in light cured composite restoration to promote the adhesion of the restoration to dental structure (dentin or enamel); In dentinal desensitization process; To bond fragments; As an adhesive cementation system complement; To the bonding of porcelain facets and repairs in porcelain.

TECHNICAL INFORMATION:

EXTRA BOND DE is a high technology bonding system with nanoparticles where primer and adhesive are supplied one-bottle (single-component). It renders an easier handling of the professional. It's a light cured