# Genius

ORIGINAL BAUMER

Indicação geral para o uso de biomateriais particulados  
General indication for the use of particulate biomaterials  
Indicación general para el uso de biomateriales particulados

BAUMER

# Genius

ORIGINAL BAUMER

## Português

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## English

Genius has a wide range of bone grafts and membranes of the market and, for more than 15 years, producing and marketing innovative, practical and safe options of biomaterials with excellent cost-benefit. Our line is composed by biomaterials of bovine and synthetic origin. All our biomaterials are biocompatible, lyophilized and sterilized, with shelf life of two years.

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<b>Indications</b>	Bone cavity.
<b>Type of Biomaterial</b>	Autogenous bone, inorganic biomaterial, biphasic biomaterial, composite biomaterial, integral biomaterial (or association of autogenous bone with some of the related materials) and membrane.
<b>Benefits</b>	Improve the contour of the alveolar process to optimize aesthetics. When only biomaterials are used there may be small bone formation. In these cases, one should use membrane as well as evaluate the periodontal biotype. When the autogenous bone associated with the biomaterial is used, there may be a layer of bone tissue.

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GenOx Inorg, GenPhos e Orthogen.

Membrana Obrigatória | Required membrane | Membrana obligatoria

GenDerm, GenDermFlex, Pratix, Techgraft.

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Carvalho PSP, Carvalho MC, Ponzoni D. Reconstruction of alveolar bone defect with autogenous bone particles and osseointegrated implants: Histologic analysis and 10 years monitoring. AnnMaxillofac Surg 2015;5:135-9.

Ho-Yeol Jang, Hyoun-Chull Kim, Sang-Chull Lee, Jang-Yeol Lee. Choice of Graft Material in Relation to Maxillary Sinus Width in Internal Sinus Floor Augmentation. J Oral Maxillofac Surg 68:1859-1868, 2010.

Jensen T, Schou S, Gundersen HJG, Forman JL, Terheyden H, Holmstrup P. Bone-to-implant contact after maxillary sinus floor augmentation with Bio-Oss and autogenous bone in different ratios in mini pigs. Clin. Oral Impl. Res. 24, 2013 / 635-644.

Thoma DS, Naenni N, Figuero E, Hämmmerle CHF, Schwarz F, Jung RE, Sanz-Sánchez I. Effects of soft tissue augmentation procedures on peri-implant health or disease: A systematic review and meta-analysis. Clin Oral Implants Res. 2018 Mar;29 Suppl 15:32-49. doi: 10.1111/cir.13114.

Wood RA, Mealey BL. Histologic Comparison of Healing After Tooth Extraction With Ridge Preservation Using Mineralized Versus Demineralized Freeze-Dried Bone Allograft. J Periodontol 2012;83:329-336.

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Prof. Paulo Sérgio Perri de Carvalho  
Prof. Alcindo Pereira Gomes Neto

## Important complementary remarks:

\* indication "off label"

\*\* biomaterials can not be exposed to the oral environment. In case it is not closed by first intention, it is recommended the use of mucoperiosteal flaps with relief in the periosteum or using the biological plug. When membranes are used, they should be buried.

\*\*\* Particulate biomaterial should not be used in areas considered contaminated. Contaminated areas are the presence of fractured roots or periodontitis. In these cases, it is recommended the extraction of the involved teeth and wait for the stabilization of the bone tissue for later surgical procedure. Surgical procedures that include horizontal and vertical enlargement may be performed with autogenous bone or with the association of autogenous bone and particulate biomaterial (50-50) in addition to the protection with fixed or screw stabilized membranes or membrane with titanium reinforcement also fixed or with titanium screen. On the titanium screen it is recommended the use of bovine pericardium membrane associated or not to PRF.

\*\*\*\* The use of a connective tissue graft associated to the filling of the gap with biomaterial is important for the maintenance of the vestibular volume, as well as in the cases of fine gingival biotype and/or a small range of keratinized tissue, it improves the prognosis long-term treatment of dental implants.

## BIBLIOGRAPHIC REFERENCES

Avila G, Hom-Lay Wang HL, Galindo-Moreno P, Misch CE, Bagramian RA, Rudek I, Benavides E, Moreno-Riestra I, Braun T, Neiva R\*. The Influence of the Bucco-Palatal Distance on Sinus Augmentation Outcomes. J Periodontol 2010;81:1041-1050.

Carvalho PSP, Rosa AL, Bassi APF, Pereira LAVD. Biomateriais aplicados a Implantodontia. ImplantNews 2010;7 (3a-PBA):56-65.

Carvalho PSP, Carvalho MC, Ponzoni D. Reconstruction of alveolar bone defect with autogenous bone particles and osseointegrated implants: Histologic analysis and 10 years monitoring. AnnMaxillofac Surg 2015;5:135-9.

Ho-Yeol Jang, Hyoun-Chull Kim, Sang-Chull Lee, Jang-Yeol Lee. Choice of Graft Material in Relation to Maxillary Sinus Width in Internal Sinus Floor Augmentation. J Oral Maxillofac Surg 68:1859-1868, 2010.

Jensen T, Schou S, Gundersen HJG, Forman JL, Terheyden H, Holmstrup P. Bone-to-implant contact after maxillary sinus floor augmentation with Bio-Oss and autogenous bone in different ratios in mini pigs. Clin. Oral Impl. Res. 24, 2013 / 635-644.

Thoma DS, Naenni N, Figuero E, Häggerle CHF, Schwarz F, Jung RE, Sanz-Sánchez I. Effects of soft tissue augmentation procedures on peri-implant health or disease: A systematic review and meta-analysis. Clin Oral Implants Res. 2018 Mar;29 Suppl 15:32-49. doi: 10.1111/cir.13114.

Wood RA, Mealey BL. Histologic Comparison of Healing After Tooth Extraction With Ridge Preservation Using Mineralized Versus Demineralized Freeze-Dried Bone Allograft. J Periodontol 2012;83:329-336.

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