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The Power Hook is a crimpable accessory designed to be attached to the working orthodontic archwire. It is manufactured from stainless steel and features a design consisting of a vertical post structure with a total length of 6.95 mm. The rectangular open base (slot) has a depth of 1.40 mm and a front height of 0.67 mm, allowing proper adaptation and fixation onto rectangular working archwires. A perpendicular stem extends from the center of the base and ends in a hook-shaped tip, designed to securely hold elastic chains, ligature wires, or pre-calibrated nickel-titanium (NiTi) coil spring rings. The distance from the center of the upper wall of the base slot to the uppermost tip of the hook measures 6.95 mm. The accessory has a smooth, polished surface with rounded edges, ensuring comfort for the patient as well as quality, safety, and reliability for the orthodontic professional.

The Power Hook is intended for use with fixed orthodontic appliances, being mechanically crimped onto the archwire at positions determined by the clinician according to the biomechanics and treatment plan. Its specific function is to serve as an anchorage or attachment point for elastic chains, intraoral intermaxillary elastics, metallic ligature wires, and calibrated nickel-titanium springs.

It facilitates space closure mechanics through anterior retraction, mesial movement of posterior teeth, or a combination of both. The Power Hook may also assist in midline correction between upper and lower arches, and in anterior intrusion mechanics. Another indication for use of the Power Hook is as a force application point for en-masse distalization of upper teeth in Class II malocclusion correction, or en-masse distalization of lower teeth in Class III correction, being associated respectively with mini-implants placed in the infrazygomatic crest or external oblique ridge regions. It may also be employed in other applications according to the clinician's treatment plan.

In mechanics involving skeletal anchorage with mini-implants, the Power Hook serves as a retention point for elastics or pre-calibrated nickel-titanium springs. Owing to its 6.95 mm hook height, it allows the application of forces with specific consideration of the center of resistance of the dental unit being moved, enabling controlled and safe tooth movements. The Power Hook is manufactured separately for the right and left sides of the dental arch.

### **Intended use**

The Power Hook Crimpable is used with fixed appliances in positions previously determined by the professional according to the mechanics to be performed and according to treatment planning, acting as a retention point or application of force through the use of a chain elastic, intraoral elastic, tie metal and pre-calibrated nickel-titanium springs, allowing the application of force, when it is desired to perform the orthodontic mechanics of closing spaces through retraction of anterior teeth, mesialization of posterior teeth, intrusion of anterior teeth, correction of midlines, mass distalization of upper teeth and mass distalization of lower teeth, respectively, in the correction of Class II and Class III malocclusions. The applicability of the crimping Power Hook is wide and can be used in other clinical situations according to professional planning.

Therefore, it has applicability in correcting malocclusions as established below:

#### **I- Sagittal problems:**

- 1) Angle Class I malocclusion;
- 2) Angle Class II malocclusion;
- 3) Angle Class III malocclusion;
- 4) Mild mandibular prognathism for compensatory treatment;
- 5) Mild mandibular retrognathism for compensatory treatment;
- 6) Smooth maxillary protrusion for compensatory treatment and
- 7) Gentle maxillary retrusion for compensatory treatment.

#### **II- Inter-arch and intra-arch problems:**

- 1) Midline deviations;
- 2) Dental positioning asymmetries.

#### **III- Vertical problems:**

- 1) Deep overbite to be corrected through intrusion of anterior teeth;
- 2) Open bites to be corrected by extrusion of anterior teeth.

### **Indication for use**

The Power Hook Crimpable is used in specific regions of the dental arch, being fixed to the working arches, together with a fixed appliance with brackets or mini-implants in positions determined by the professional according to the mechanics to be used and according to treatment planning. They serve as an anchor or retention point for a chain elastic, intraoral elastic, metallic tie, pre-calibrated nickel-titanium spring, to carry out the mechanics mentioned above and others in which they may have applications.

Clinical indications - Deviations from ideal occlusion (mentioned above in Intended Use) and the associated psychosocial issues today play an important role in defining the need for orthodontic treatment.

### **How to use:**

The Power Hook Crimpable is fixed to the rectangular arch by means of mechanical constriction of the fixation base using crimping pliers, Weingart or similar in positions determined by biomechanics and treatment objectives. To avoid the risk of the accessory sliding due to the application of force, the attachment point can be reinforced using a spot welding machine with a special electrode, or even using the welding feature with silver solder and flux. Correct use of the Power Hook Crimpable involves the following steps:

- 1) Selection of the appropriate orthodontic arch;
- 2) Adaptation of the Power Hook Crimpable to the rectangular arch in the position predetermined by the mechanics to be used;
- 3) Press the accessory fixing base using crimping pliers, Weingart or similar;
- 4) The attachment of the accessory to the arch is reinforced with spot welding or silver solder and a torch, these two resources being optional;
- 5) Finishing and polishing the welding area.

### **Removal of the Power Hook Crimpable:**

After achieving the proposed treatment objectives, the accessory must be removed and the following steps must be followed:

- 1) Remove the arch from the oral cavity;
- 2) Using 139 or Weingart pliers positioned perpendicular to the fixing base, apply pressure;
- 3) The Power Hook Crimpable can be removed.

The professional may also choose to eliminate the working arch with the orthodontic accessory and decide to create a new arch.

The Power Hook Crimpable may undergo changes in their surface structure when exposed to adverse conditions present in the oral environment. In the oral cavity, the biodegradation of metals usually occurs through electrochemical breakdown. Although electrolytic attack appears to be the major cause of the derangement, bacteria and their cellular degradation products, and selective interactions with gases such as oxygen and carbon dioxide may contribute to the degradation of materials in the oral cavity. Organic acids and enzymes in particular can also affect the various metals and the pH of the environment where the devices are installed, exerting a great influence on the corrosive process.

The medical product in question can only be placed on a patient by a dental surgeon specialist in Orthodontics.

**Composition:**

Power Hook Crimpable – AISI 304 stainless steel

**Storage Conditions:**

In intact and closed packaging.

This product does not require temperature, humidity, or light control during storage. The recommendation is to keep it in a clean, dry place.

The manufacturer recommends the single use of this product. It should be disposed of correctly in accordance with the information contained in the Manual of Dental Services Risk Prevention and Control, Technology in Health Services, 2006 (available on the ANVISA website).

**Warnings and Symbology:**

- It is essential that only the specialist installs the product on the patient, since in addition to the clinical examination, the professional must have the necessary knowledge to make a systematic description of the patient's dentofacial features, list the priority problems, describe possible solutions for each of the problems, considering the possible favorable or unfavorable interactions for each of the proposed attempts.
- It is also the orthodontist's duty to determine the maximum expected gain from orthodontic therapy and the minimum loss, or to calculate the benefit versus risk ratio for each patient, with both risk and benefit including psychosomatic factors to determine the best management for each patient.
- Patients must be instructed about oral hygiene, maintenance and care with food. The following should be avoided: ice, coconut, hard and/or sticky sweets, caramels, nougats, jawbreakers, corn, popcorn and sugar cane, chewing gum and caramels. Biting hard objects such as pens and pencils produces the same effects as the foods mentioned above and can displace the Power Hook Crimpable from its fixing position or even bend the vertical rod. Therefore, patients should be advised to avoid hard foods.
- This product should not be used in an MRI environment, as it may cause interference in the images during the examination. It is up to the orthodontic professional to warn the patient of this contraindication and arrange for the appliance to be removed if necessary.
- The high concentration of bacteria and fungi present in the mouth can accelerate the corrosion of metal orthodontic appliances. The pH of the oral environment also has a major influence on the corrosive process, so a diet rich in acidic foods, such as excessive intake of soft drinks with acidic pHs, should be avoided as much as possible during treatment.
- Do not touch the Power Hook Crimpable with ultrasonic tips: this can lead to fractures and detachment.
- The product is supplied in non-sterile condition, however it has proven biocompatibility and no cleaning, disinfection or sterilization measures are required, provided the packaging has not been tampered with. If the product comes into contact with surfaces or substances other than those intended for use, it is recommended that it be discarded.
- The product must not be reused or reprocessed. Do not use the product if its packaging is damaged

or open, to minimize possible dangers related to the transmission of microorganisms or possible product counterfeiting.

- It is the responsibility of qualified professionals to advise their patients of the risk of aspiration and swallowing in the event of breakage/ detachment of parts, not to underestimate the possible immediate or delayed complications resulting from this type of accident and to adequately inform the patient and/or their legal guardians of the risks inherent in the event and the measures to be taken.
- This product should only be used by qualified orthodontic professionals. Use by people without the necessary technical knowledge could lead to unwanted tooth movement, bone loss, tooth loss, fenestration, gingival recession, or root dehiscence.
- It is the responsibility of the qualified professional to advise the patient on correct oral hygiene, in order to avoid the appearance of plaque and tartar or diseases such as gingival inflammation (gingivitis), periodontitis or even endocarditis. Considering that the product is used in contact with mucous membranes and body fluids, it is recommended that professionals apply good clinical practices for the disposal of medical products in accordance with current health regulations.

**ATTENTION: Orthometric cannot be held responsible for the results obtained by applying the product without the necessary precautions or failing to observe the warnings.**

### Contraindications

Orthometric only supplies its products to qualified professionals. It is the orthodontist's sole responsibility to identify conditions that may be contraindicated for treatment, such as:

- Patients with poor oral hygiene.
- Patients unable to cooperate with treatment.
- Patients with an oral environment conducive to the demineralization of tooth enamel.

	Single-use product. Do not re-use		Caution
	Consult instructions for use		Non-sterile
	Medical device		Contains Nickel
	Catalogue number		Date of manufacture

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Revisão / Revision: 0  
Data / Date: 13/11/2025

	Batch code		Keep away from sunlight
	Manufacturer		Do not use if package is damaged
	Use-by date (if packaging is unopened)		Keep dry on during transport and storage.