14089 Sigma Flow - Instr de Uso PEI

ESPAÑOL

SIGMA FLOW Composite Fluido

EXCLUSIVO PARA USO PROFESIONAL

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510mm

ENGLISH

SIGMA FLOW Flowable Composite Resin

PROFESSIONAL USE ONLY

SIGMA FLOW is a composite resin of low viscosity, which provides a good resistance to compression, diametrical tension and bending.

Composition:

SIGMA FLOW contains: bis-GMA, Dimethacrylate Resins, Glass Boron-Aluminum Silicate, Synthetic Silica and Dyes.

Presentations:

Box containing: 1 syringe with 2.4g of SIGMAFLOW + 1 application tip.

Indications:

SIGMA FLOW is indicated for:

- Dental restoration for enamel hypoplasia, fluorosis, erosion and abrasion injuries or pigmented.
- Restorations Classes III and V.
- Sealing of pits and fissures.
- Repair of composite restorations and porcelain.
- Filling of excavated areas in class II restorations.
- Cementation of porcelain veneers.
- Tooth fragment bonding.
- Base or liner for direct and indirect restorations.

Contraindications:

SIGMA FLOW é is contraindicated in patients with a history of sensitivity to any component of the formulation or composites with acrylates. If allergic reaction occurs, seek medical attention, remove the product if necessary and discontinue its use.

Warnings:

- The use of gloves is recommended during the utilization of SIGMA FLOW.
- Avoid using SIGMA FLOW in patients with a history of allergy to composite resins with acrylates. Do not use the product in patients with a history of allergy to the material. In case of allergy, immediately discontinue use and advise the patient to seek medical advice. Avoid the product contact with oral tissue or skin. In case of contact, immediately rinse the area with water. Avoid contact with the eyes; in case of accidental contact, wash the eyes immediately with water and see a

doctor. This product must be used only by dentists and for the intended use described in the Instructions for Use.

- Make sure the application tip is unobstructed and the syringe working well before starting the application of the product.
- Avoid using eugenol-based materials because they prevent the polymerization of SIGMA FLOW.
- After opening the product it should be stored in its original package tightly closed, under the conditions described in the item "Storage and conservation", otherwise the product will lose its characteristics.
- Wastes or residues of the used material lose their characteristics and therefore should not be reused. They must be disposed as the local law in force.
- Keep out of the reach of children.
- SIGMA FLOW is intended only for professional use.

Instructions for Use:

"Attention when using a new applicator tip: Gently squeeze the plunger and wait for the resin extrusion. If you apply excessive force, the material extrusion will continue even after the digital pressure is finished."

For restorations:

- 1- Select color with hydrated teeth and under natural light.
- 2- Perform preparation for adhesive restorations.
- 3- If necessary, use calcium hydroxide for pulp protection.
- 4- Make insulation and apply the acid conditioner ALPHA ETCH on enamel for 30 seconds and dentin for 15 seconds. Rinse and remove excess of moisture avoiding dehydrating the dentin.
- 5- Apply SIGMA FLOW with the aid of the applicator tip. Use a brush to adapt the material and photoactivate for 30 seconds. In the case of severely pigmented enamel, first apply a layer of opaque SIGMA FLOW OPACO and photoactivate for 30 seconds.

NOTE: In case of Class V restorations or little retentive cavities, the use of a dentine adhesive (e.g. ALFA BOND DE PLUS) is recommended to increase the adhesiveness of the resin.

- 6- Follow the conventional finishing and polishing recommended procedures for light-cured resins.

For filling of excavated areas:

*Follow steps 1 to 5 described above.

- 6- Apply SIGMA FLOW to fill and level the excavated areas of difficult access. Do not allow the product extends the occlusal region. Photoactivate for 30 seconds.
- 7- Proceed insertion of conventional resin in areas easily accessible and of occlusal contact.

- 8- Follow the conventional finishing and polishing recommended procedures for light-cured resins.

For Sealing of pits and fissures:

- 1- Perform the absolute isolation and acid etching of the surface of pits and fissures. Rinse and dry.
- 2- Apply SIGMA FLOW with the aid of the applicator tip. Do not allow the product extends the occlusal region. Photoactivate for 30 seconds.
- 3- Follow the conventional polishing recommended procedures for light-cured resins.

For Cementation of porcelain veneers:

- 1- Select the right color by a test of veneer with SIGMA FLOW in the tooth itself, without photoactivate.
- 2- Perform the internal conditioning of ceramic with hydrofluoric acid. Rinse and dry. Apply silane in the conditioned surface.
- 3- perform isolation and acid conditioning of the enamel . Apply the adhesive ALFA BOND DE PLUS without photoactivate .
- 4- Apply SIGMA FLOW on the inner surface of the workpiece and place the veneer. Remove excess.
- 5- Photoactivate for 30 seconds.

Storage and Conservation:

SIGMA FLOW must be stored between 2°C and 30°C in its original package. Keep the bottles tightly closed after use. For better conservation keep it in refrigerator. Do not expose to excessive heat. Under these conditions, the product is valid for up to 2 years as of the manufacturing date. Do not use this product if shelf life had expired.

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DFL Indústria e Comércio S.A.