CONTEÚDO:

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CONTENT:

1 vial containing 520 mg
Composition: Poly-L-lactic Acid, Carboxymethyl Cellulose (CMC) and Mannitol
Weight after reconstitution: 16.52 g
Volume after reconstitution: 16 ml
Each vial with lyophilized powder contains:
Poly-L-lactic Acid 210 mg
Carboxymethylcellulose 132 mg
Mannitol 178 mg

FEATURES

Provided as a lyophilized sterile preparation for injection in amber glass vial that is sealed by a penetrable lid, covered by aluminum seal with removable lid.
• Package unit: 1 vial/box
• Single use
• Shelf life: 2 years

INTENDED USE DESCRIPTION

Rennaova Elleva® consists of Poly-L-lactic Acid in sterile suspension, which is reconstituted from sterile lyophilized powder adding sterilized water for injection.
Rennaova Elleva® should be injected in deep dermis, subcutaneous or supraparietal area to improve skin flaccidity arisen from intrinsic and extrinsic aging process and volumetric correction of depressed areas, as grooves, wrinkles, skin depressions and atrophic scars in the face, neck, cleavage, arms, abdomen, thighs and buttocks. Facial changes in regions such as, temporal, malar and zygomatic, due to fat atrophy (lipoatrophy), mainly from facial fat absorption in patients HIV positive, can also be remedied by Rennaova Elleva® injection, in deep dermis, subcutaneous or supraparietal area. Rennaova Elleva® is gradually absorbed by the body naturally and fully absorbed after 18 months.

VOLUME PER INJECTION

Rennaova Elleva® maximum volume for each individual injection should be limited to 0.14 ml – 0.28 ml, with a 0.5 to 1 cm distance apart. Avoid overcorrection. Recommendation is to apply up to 0.7 ml every cm², 3 - 6 session, with 22.8 ml product total volume.

PREPARATION BEFORE USE

1. Remove the vial tear-off cap and clean the vial penetrable stopper with antiseptic. If the vial, seal or tear-off cap are damaged, do not use.
2. Attach 18G sterile needle (recommended) to a 20 ml sterile use syringe.
3. Introduce the sterile needle in the vial stopper and slowly add 16 ml sterilized water (recommended) in the vial.
4. Stir the vial for 10 minutes. Allow the vial to stand at least for 1 hour to ensure full hydration, attaining uniform translucent suspension. Immediately stir again before use. A single vial stirrer can be used.
5. Reconstituted product is for single and individual use, and can be used between 1 to 72 hours after reconstitution, stored under 2-8°C temperature.
6. Dispose any remaining material after use or 72 hours after reconstitution.
7. Treatment areas should be periodically massaged to distribute the product evenly.
8. Introduce a straight, sterile and bevel-up needle in the skin with approximately 30 to 40 degrees angle to face areas and 45 to 60 degrees in the other areas related to the skin, and, following, insert needle till reaching the skin desired deepness. Change in tissue resistance is noted when the needle crosses dermis to subcutaneous layer. If the needle is inserted in a very superficial angle, or if the needle tip is not enough inserted, the needle tip may be in the medium or superficial (papillary) dermis. If the product is also injected superficially, the injected area will immediately blanch or briefly after the injection. If the area injected blanches, remove the needle and massage the area. In case of blanching not disappearing, no further injection can be carried out to the patient.
9. Normally, patients undergo some edema degree related to the injection procedure itself, which will give the appearance of a full correction till the end of injection session (within 30 minutes). Patient should be informed that the edema related to the injection generally disappears after some hours or after some days, resulting in the original contour reappearance.

INSTRUCTIONS OF USE

1. Before the treatment with this product, a patient should receive patient labeling and be fully informed by the health care professional skilled about the use intended, indications of use, as well as contraindications, use warnings and precautions, correction expected and possible side effects and care after the procedure. Each patient should be informed that this product quantity and number of injection sessions will depend on the patient's need. In addition, for desired outcomes, patients should be informed that more than one session is needed.
2. Each injection session should be carried out with aseptic technique and universal precautions. The patient should be informed about the potential contact with the patient's body fluids: blood, from the injection site. A topical anesthetic can be used in the application site to reduce sensitivity and pain during the treatment.
3. Before injecting in subcutaneous layer, 16 ml sterile water for reconstitution should be added to each use. This product is to be injected only into dermis and subcutaneously layer.
4. To the injection 21G to 26G disposable needles and/or 22G to 25G disposable micro cannula can be used. To keep the suspension uniform during the procedure, intermittently stir the product in the syringe. Before initial injection, expel some drops through the 21G needle fixed to remove the air and check for needle clogging. If the needle is dogged during injection session, replacement is necessary. In case of clogging, remove the needle, remove a small quantity of the product, attach a new sterile needle, then withdraw some drops to remove the air and check for needle clogging again.
5. Treatment areas should be periodically massaged to distribute the product evenly.
6. Introduce a straight, sterile and bevel-up needle in the skin with approximately 30 to 40 degrees angle to face areas and 45 to 60 degrees in the other areas related to the skin, and, following, insert needle till reaching the skin desired deepness. Change in tissue resistance is noted when the needle crosses dermis to subcutaneous layer. If the needle is inserted in a very superficial angle, or if the needle tip is not enough inserted, the needle tip may be in the medium or superficial (papillary) dermis. If the product is also injected superficially, the injected area will immediately blanch or briefly after the injection. If the area injected blanches, remove the needle and massage the area. In case of blanching not disappearing, no further injection can be carried out to the patient.
7. Normally, patients undergo some edema degree related to the injection procedure itself, which will give the appearance of a full correction till the end of injection session (within 30 minutes). Patient should be informed that the edema related to the injection generally disappears after some hours or after some days, resulting in the original contour reappearance.
8. Immediately after this product injection session, redness, swelling and/or hematoma may occur in the treatment area. After the injection session, ice pack should be applied (avoiding any ice direct contact with the skin) in areas treated to reduce edema and/or bruise. To help evenly distribution of the product in contour deficiency, it is important to massage by hand the area treated for 5 minutes and 5 times daily, at least for 5 days. Physical activities should be performed according to medical advice. This product is visible through ultrasound image and magnetic resonance, but it is not visible through computed tomography or radiography.
9. Exam the patient again within 4 weeks, after the injection session to determine if additional correction is needed. Original contour deficiency may reappear initially, but it is expected the deficiency improves gradually within some weeks. Patient should be

informed on the additional injection sessions potential need in the first appointment.
10. During the first injection session, only limited correction should be carried out. Area treated should never be excessively corrected (sobrecorreção) just in one session.
11. This product powder form should be stored under room temperature (5–20°C) avoiding direct sunlight. Refrigeration is not necessary. After reconstitution, the product must be stored under 2-8°C temperature and can be used within 72 hours.

CONTRAINDICATIONS

1. This product should not be used in a person suffering from active inflammation or infection in or near the treatment area.
2. This product should not be used in a person hypersensitive to any of its components.
3. Do not use in periorbital area, as well as its proximities. Do not use on the lips.
4. This product should not be administered in superficial and medium dermis, as well as to thin wrinkles correction/treatment.

WARNINGS

1. This product should only be used in deep dermis or in subcutaneous layer.
2. To avoid papules and nodules near the treatment area, superficial injections should be avoided. To mitigate papules and nodules, the area treated should be massaged so the product is evenly distributed.
3. To avoid contamination, each vial is for the patient's individual use. Do not reuse nor sterilize the vial again. Do not use if the package is damaged or opened.
4. This product should not be implanted in blood vessels. This product implantation in dermal vessels may cause vascular occlusion, infarction or embolic phenomena.
5. Do not over correct contour deficiency, as the depression should gradually improve during several weeks after the injection, as the treatment effect occurs.
6. Do not inject into the red (middle) area of the lip.
7. Powder should be reconstituted with sterilized water and can be used within 72 hours when stored between 2-8°C.
8. Sterile needle and syringe must be used to inject the product.

PRECAUTIONS

1. This product should only be used by a health care professional skilled who is fully familiarized with the product and with the product education materials.
2. Treatment area should be cleaned with antiseptic and must not be inflamed or infected.
3. Such as occurs with all injections, patients with coagulation failures or simultaneously under anticoagulant therapy are under greater risk of hematoma, bruise and/or bleeding formation in the injection site.
4. This product safety and effectiveness have not been assessed in pregnant, nursing patients or less than 18 years old.
5. Interaction studies were not carried out of this product with drugs or other substances or implants.
6. High risk of papules and nodules was reported in published literature after injections in periorbital area.
7. This product safety and effectiveness were not assessed in patients with keloid or hypertrophic scars formation. This product should not be used in patients with history known or susceptibility to keloid or hypertrophic scars formation.
8. The patient should be informed that he or she should reduce the treatment area exposure to the sun and avoid UV lamp exposure, till all swelling and initial hematoma have been settled.

ADVERSE REACTION

The following adverse reactions were found as result of clinical studies, post market surveillance and/or academic papers.
• Reaction around the area treated
- Bleeding
- Pain
- Induration
- Swelling
- Hematoma
- Ecchymosis
• Immunological disorder
- Hypersensitivity
- Angioedema
- Skin sarcoidosis
• Infection and Inflammation
- Cellulitis e Staphylococcal infection
- Abscess
• Skin and subcutaneous tissue disorder
- Bruise
- Hematoma
- Skin atrophy and hypertrophy
- Erythema
- Sclerotherapy
• Majority of papules are not visible. Although they can be intensified, there are no symptoms.
- Granuloma
- Scar
- Skin blemish
- Nodules

STERILE. Sterilized by gamma radiation. Single use product.
REPROCESSING PROHIBITED.

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