## Bio-indicador 24 horas

### PORTEGUÉS

#### Conteúdo:

Caixa com 10 frascos plásticos contendo tiras de papel com esporos Geobacillus stearothermophilus ATCC 7953 e ampolas de vidro seladas.

#### Indicação:

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ANSI/AAMI ST:59:1999, Sterilization of health care products Biological Indicators

Part 1: General.

ANSI/AAMI ST:19:1999, Sterilization of health care products Biological Indicators

Part 3: Biological Indicators for moist heat sterilization.

ISO 11138-1:1994, Sterilization of health care products Biological Indicators

Part 1 General.

ANSI/AAMI ST:19:1999, Sterilization of health care products Biological Indicators

Part 3: Biological Indicators for moist heat sterilization.

ISO 11138-1:1994, Sterilization of health care products Biological Indicators

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### ENGLISH

#### Contents:

Box with 10 bottles containing paper strips with Geobacillus stearothermophilus ATCC 7953 spores and sealed glass ampoules.

#### Indication:

Indicated for periodical monitoring of sterilization cycles and in validation of autoclaves.

#### Product Description:

Sterilization is essential in the control of infections in the various health establishments. The BIOINDICATOR aims to monitor steam sterilization, providing greater reliability in this process, through bacterial spores resistant to heat, the Geobacillus stearothermophilus ATCC 7953.

#### Composition:

Culture medium (Tryptic soy broth) and pH indicator-Bromocresol purple; Geobacillus stearothermophilus ATCC 7953: minimum population of 105 to 106 spores inoculated in paper disk.

#### Mode of use:

1. Remove 02 ampoules from the box and identify in each label the equipment, charge and date. Place only 01 ampoule in the paper envelope surgical grade. The other ampoule stay out of the autoclave, it is the control ampoule.
2. Place the paper envelope surgical grade in the autoclave in a local of difficult access, for example in a lower part/bottom.
3. Start the cycle according the instruction of use of the manufacturer. After finalized remove the envelope from the autoclave and wait it cool down for 10 minutes.
4. Remove the Bio Indicator of the envelope. To break it compress the tube against the wall of the incubator until the internal glass ampoule breaks and its culture medium is released to come into contact with the spores that are in the paper disc. Tap slightly with the index finger. Be careful to not brake the plastic tube and keep the tube always in vertical position. It can also use the ampoule breaker.
5. Repeat this same operation of breaking with the control ampoule.
6. Start the process of incubation of 24 hours with temperature of 55 to 60°C.
7. Make the results registers of the two ampoules in the controlling sheets, according the model shown.

#### Results analysis:

**Test ampoule:** if the changing of color doesn't occur (negative result) indicates that there's no bacterial growth, in other words, it has occurred the destruction of the spores and the medium remain purple. This tube was sterilized correctly. If occurs the color changing to yellow positive result, in other words, it has occurred bacterial growth indicating that the sterilization was inappropriate. Interrupt the use of the autoclave, analyze the physical parameters, if there is overloading of the equipment and repeat the cycle. If the color yellow remain (positive results), forward the equipment to technical assistance. Upon returning, repeat the procedure.

**Control ampoule:** this tube does not undergo any sterilization process, so their turning yellow (positive results) should occur generally before 4 hours of incubation. After turning point, register the result, autoclave the tube and discard. If after the incubation of this tube the result was negative (without color changing) this result invalidate the equipment test and it can occur due malfunction of the incubator, inadequate tube sterilization, inadequate breaking of the glass ampoule or inappropriate storage of the product.

**Results:** the test is considered valid or approved if the two ampoules were incubated together and when the test ampoule obtain negative results and the control ampoule positive result. Any result different than that the test should be reproved.

#### Precautions:

- Strictly follow the procedure provided in the mode of use;
- Start the process of incubation no more than two hours after sterilization;
- In case it is used more than one test ampoule, it can be used only one control ampoule;
- The control ampoules should be with the same lot and fabrication of the test ampoule;
- Check the label to the ampoule if occurred any changing in the chemical indicator from pink to brown, after sterilization;
- Do not reuse the ampoules;
- Do the discard of the ampoules only after sterilized in autoclaves;
- The Bio Indicator does not withstand temperature exceeding 130°C, and so it must not be used in greenhouses or sterilization ovens.

#### Storage and Transportation Conditions:

Store out of the sun, at a temperature between 15 and 30° C and relative humidity between 35 and 60%.

#### Mode of disposal:

All the material used in the microbiological control must be autoclaved for at least 30 minutes at a temperature of 121° C, and then discarded in hospital trash.

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ANSI/AAMI ST:59:1999, **Sterilization of health care products Biological Indicators**

Part 1: General.

ANSI/AAMI ST:19:1999, **Sterilization of health care products Biological Indicators**

Part 3: Biological Indicators for moist heat sterilization.

ISO 11138-1:1994, **Sterilization of health care products Biological Indicators**

Part 1 General.

#### Biological indicators of resistance data Bio-indicator for sterilization Steam:

Population (nº average spores / unit), survival time in minutes, time of death in minutes and Amount D (121°C á 134° C) it is determined with each batch. Information they are on the adhesive label.

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Odontológicos S.A.**