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
Microcut

Microcut is a hack-saw frame to which saws or interchangeable diamond strips are adapted. It is a device designed to remove excesses of composites, adhesives, cements, and amalgams from interproximal areas without damaging the contact points.

Use: The Microcut hack-saw frame move longitudinally, in back and forth movements within the inter-proximal space applying low force against the excess to be removed. With the diamond strip installed in the hack-saw frame, the reduction of the existing excesses is accomplished, as well as the final finishing in the approximal area.

Blade replacement: loosen partially the nut counterclockwise (1) and remove the used blade. Insert the new blade (2), firstly fixing one end to the arc hitch, and making its contour (3) then, hitch the other end of the blade, following the same fixation process. Retighten the nut (4) clockwise. Make sure that the nut is firmly tightened so that the system (arc + saw or arc + strip) works properly. Attention: do not make wrinkles on the blades to fit the arc. The blades, due to its thickness, can break with this movement.

Warnings, precautions and restrictions: Professional use only. The saw blade and diamond strip are single use devices. Its reuse can cause cross contamination. The saw blade edges are extremely aggressive. We recommend the protection of the papilla, for example, applying a wedge to avoid accidents. The hack-saw frame should be completely disassembled for a proper hygienization and sterilization. Before sterilization, proceed with the cleaning using water and enzymatic detergent. When there is proven allergy to any of its components, the product must not be used. To avoid an accidental swallowing of this product it is necessary the application of a total isolation technique of the operatory field. The product is a medical device. Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Sterilization: The hack-saw frame may be sterilized by autoclave (126°C - 16min). Before sterilization, proceed with the cleaning using water and enzymatic detergent.

Composition: Hack-saw Frame + saw blade: Stainless steel. Diamond Strip: Stainless steel and diamond.

Storage: To be stored in a dry place, sheltered from sunlight.

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