ENGLISH

Presentation: 1 syringe with 4 g end 1 instruction manual.

Description: FIX BITE is indicated for temporary malocclusion of the teeth during the orthodontic treatment. Due to its high adhesiveness to the tooth enamel, the product is also indicated for cementing orthodontic bands.

Indication: The product is applied in order to apply bite raisings which will "unlock" the occlusal contact and will allow the transverse correction of the patient's bite.

Directions of use:

For bite raising: 1. Select the tooth and perform the proper prophylaxis of its occlusal surface. After having cleaned it, dry it properly; **2.** By using 37% phosphoric acid, perform the conditioning of the region for 37% seconds. Rinse, dry and isolate the region properly through an absolute isolation or an efficient relative isolation; **3.** Apply Bite on increments of 2 mm maximum and photopolymerize them by 20 seconds with the help of a photopolymerizer with a minimum power of 400 mW/cm². Add increments until the desired raising is reached; **4.** Check the level of the occlusal raising with the help of an articulating paper. Generally, the level of the raising is the lowest as possible in order not to cause a malocclusion of about 1 mm on the other teeth. The excess of resin may lead to excessively strong occlusal contacts which may cause occlusal dysfunctions and tooth pain to the patient; **5.** Since it is a product for temporary use, its removal must be immediately performed after the planned period for treatment. Use multi laminated drills. **Note 1:** Additional increments to the raising may be added, if necessary, by following the procedure below: **a)** Dry the area and perform the relative isolation;

b) Apply the first coat of adhesive in the region which will receive the additional increments and photopolymerize it by 10 seconds, apply a second coat and photopolymerize it by 10 seconds. Apply smoothly a compressed air jet by 10 seconds and then polymerize by 10 seconds. **c)** Add the increments as required, photopolymerizing them by 20 seconds. **Note 2:** The application of Bite does not require the use of adhesive when the enamel is conditioned with 37% phosphoric acid. In this condition, the performance of this material is effective for most patients. However, in cases which the resin gets an excessive occlusal effort or when the tooth presents an occlusal surface with deep, irregular grooves, the use of adhesive, applied before the resin, is recommended, due to the increase of resistance this product offers. Apply the adhesive according to the instructions of your supplier.

For cementation of bands: 1. Perform the prophylaxis of the tooth which will receive the orthodontic band; **2.** By using 37% phosphoric acid, condition the free surfaces of the teeth (vestibular and lingual/palatal) by 30 seconds. Rinse and dry the surface; **3.** Use a relative isolation to properly isolate the area; **4.** Asperse the interior of the orthodontic band to increase the adhesion; **5.** Apply Bite on the interior of the band and perform its positioning on the tooth, make sure the positioning is correct; **6.** Use a dental explorer or a dental spatula to remove the excess of the product (mainly on occlusal and cervical surfaces); **7.** Use a photopolymerizer (400 mW/cm² minimum power) to photoactivate the occlusal surface by 40 seconds and the cervical surface by 40 seconds. **Note:** The patient must abstain to food by, at least, 30 minutes after the cementation of the orthodontic band.

Composition:

Uretano Dimetacrilato
Trietilenoglicol Dimetacrilato
Fosfato de Hidroxietil Metacrilato
Canforquinona
Ethyl p-(dimethylamino)benzoate
Butylated hydroxytoluene (BHT)
hydroxyphenol benzotriazole
Dióxido de silício coloidal

Storage conditions: store the syringe in a fresh, dry place, protected from light. keep in a temperature range from 5° to 30° c. do not keep near materials, which contain eugenol or clove oil in its composition.

Warning: Exclusive use of the dental surgeon in case of allergic reactions, suspend the use of the product. Products such as eugenol or clove oil inhibit the polymerization. Close the syringe immediately after use. The use of cationic solutions, plaque disclosing and chlorhexidine may cause discoloration of the resin. Do not use the product if it is expired.

Precautions: Close the syringe immediately after use. Use goggles while handling the photopolymerizer. Keep the product away from children. In case of accidental contact with the eyes, rinse the area with water and call a doctor immediately.

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