## FIBER POST

## FIBER POST

# 0,5

20,0 x 1,4 x 0,65mm

# 0,5

20,0 x 1,4 x 0,65mm

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**DESCRIPTION:** Fiber Post is an intraradicular reinforcing post manufactured with a matrix based on fiberglass and epoxy resin of high mechanical strength. Fiber Post is esthetic, translucent and radiopaque. Its dual taper allows better adaptation to root canals without demanding greater wear of the structure. It can be used in direct or indirect restorations in anterior and posterior teeth.

**COMPOSITION:** Fiberglass and epoxy resin.

**INDICATION:** Indicated for reconstruction of teeth with extensive coronary destruction and satisfactory endodontic treatment. It aids in the support and anchoring of restorations, providing retention and resistance to definitive restorative material (cases of direct restorations) or in core build-up (for indirect restorations).

### POST SIZE:

<b>0,5/1</b>	<ol style="list-style-type: none"> <li>Molars and pre molars;</li> <li>Central and lateral incisors;</li> <li>Case of partial loss of tooth structure;</li> <li>Reinforcement of nonvital teeth, weakened by excessive canal wear.</li> </ol>
<b>2/3</b>	<ol style="list-style-type: none"> <li>Cases with greater mechanical stress (significant loss of tooth structure);</li> <li>Larger canals (with higher wear);</li> <li>Reinforcement of nonvital teeth, weakened by excessive canal wear.</li> </ol>

## ENGLISH

**HOW TO USE:** 1. Perform the radiographic examination (periapical radiography) and check the quality of endodontic treatment and periapical health; 2. Select the post according to the conduit diameter by positioning it on the radiograph of the tooth; 3. Perform absolute isolation of the operative field; 4. Remove endodontic filling of 2/3 the root canal, using suitable instruments. Leave at least 3 to 5mm of gutta-percha in the apical region, sealing and ensuring the prevention of bacterial contamination; 5. Prepare the root canal with Largo drills corresponding to the diameter of the post. Avoid unnecessary wear; 6. Make a periapical radiograph of the selected post within the root canal. Check in the image preparation of the conduit, the post adjustment and the condition of the remaining filling. The post must be perfectly adapted to the executed preparation, with no spaces between the end of the post and the filling of the apical remnant; 7. If necessary, cut the coronal surplus post with diamond drill in high rotation under refrigeration. Return to position for adjustment check; 8. Clean the post with 70% alcohol and dry; 9. Apply a layer of Silano Maquira on the entire surface of the post and allow to dry for 1 minute. The post should be completely dry at the time of cementation; 10. Apply Phosphoric Acid Gel 37% Maquira for 15 seconds inside the root canal and remaining coronary. Thoroughly wash and remove excess water from the canal with absorbent paper cones; 11. Apply the adhesive system (according to the manufacturer's instructions) across the full extent of the conditioned root canal and allow to photopolymerize; 12. Carry the resin cement in the root canal using appropriate instruments. If you prefer, cover the post with the cementing agent too; 13. Insert the post into the canal. Note that the post must be inserted to the bottom of the prepared cavity and the intracanal portion of the post should be completely covered by the resin cement; 14. Remove excess cement; 15. Carry out the photopolymerization of the cement (according to manufacturer's specification) through the surface and through the post; 16. It is recommended to perform an x-ray to evaluate the final quality of the post cementation.

### PREPARATION SEQUENCE WITH LARGO DRILL:

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### CONTRAINDICATIONS:

Fiber Post is contraindicated to the following patients:

- Patient with parafunction.
- Patient that has a coronary remnant less than 2 mm.
- Patient with dental element not treated endodontically or with unsatisfactory pulp therapy.
- Patient with hypersensitivity to the components of the product.

**PRECAUTIONS AND WARNINGS:** 1. This product is intended for dental use only and must be handled in accordance with its instructions for use. The manufacturer is not liable for any damage caused by off label use or incorrect manipulation. 2. Restricted use for professionals specialized in dentistry, according to local legislation and with specific training for their manipulation. 3. The indication for use of the endodontic post depends on careful evaluation of the quality of the endodontic treatment, periodontal support and periapical health of the dental element. 4. Proper selection of post size is critical to expected clinical results. 5. Do not use the product if it is out of date. For product disposal follow the applicable local, state and/or federal regulations. 6. Handling of the product must be carried out under absolute isolation. 7. Fiberglass posts may present exposed tips of fiberglass which can puncture the skin. Cutting the posts to adjust the length can also expose sharp fibers.

## ENGLISH

## ENGLISH

- Use PPE when handling the product, including operator gloves and safety glasses, as fragments of the product may come loose during this process.
- Wear a mask and avoid inhaling the powder from the cut.
- Use forceps to hold the post.
- The cutting of the fiberglass post must be done with diamond tips in high rotation and under refrigeration. Do not use scissors or pliers.
- Avoid cutting it close to the patient's face to avoid any accidents or discomfort with cutting residues.
- Avoid breathing the dust from cutting the posts.
- The post must be free of oily soils, organic materials or any other contaminating agents at the time of cementation to avoid adhesion failure.
- The use of dentin adhesives over Fiber Post is unnecessary and not recommended, because the adhesive is a polymer with mechanical properties that are inferior to those of the post and resin cement, which produces an interface of low adhesive strength between the components. During the process of cementing the post previously silanized with Silano Maquira, direct contact with the cement promotes a gradient of favorable mechanical properties.
- The use of fiberglass post is not indicated in rehabilitations of patients with parafunctional habits.
- The technique of core build-up with post in fiberglass should only be used on teeth with a minimum of 2 mm of height of coronary dentin.
- Complete removal of the endodontic filling and absence of apical seal may cause bacterial contamination of the root canal.
- The fiber post must be completely covered by the cementing agent or by the composite resin.
- The product should not be used in patients with hypersensitivity to the components of the product.
- In case of accidental swallowing or aspiration of the product, seek medical attention immediately.
- Keep the product at room temperature and its original packaging properly closed.
- Keep the product out of the reach of children.
- Single use product.
- Prohibited reprocessing. Reprocessed may cause change in product properties, post failure, loss of restoration, contamination or root fracture.

**CONSERVATION AND STORAGE:** Keep the product at room temperature in its original packaging always well sealed.