# REFORPOST®

## PORTEGUÉS

### INDICAÇÕES

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## ENGLISH

### INDICATIONS

For the support of restorations and prosthetic crowns.

### COMPOSITION

Peeso drills.....	Stainless steel
Rubber Stops.....	Latex
Template for radiographic measurement.....	Plastic
<b>REFORPOST® GLASS FIBER:</b>	Glass fiber (80%); pigmented resin (19%), stainless steel filament (1%).
<b>REFORPOST® CARBON FIBER:</b>	Carbon fiber (72%); pigmented resin (22%), stainless steel filament (6%).

### TECHNIQUE OF USE

#### A. ROOT CANAL PREPARATION

1. Select the pin using the template and the periapical x-ray of the tooth (photo 1);
2. Remove some of the root canal filling with heated instruments to the predetermined depth. The use of drills can cause apical displacement of the gutta-percha. Allow at least 3 to 5 mm of filling at the root apex due to the presence of accessory canals, which must remain sealed to prevent bacterial contamination (photo 2);
3. Prepare the canal with the drill corresponding to the diameter of the selected post (photo 3).

#### B. PREPARATION OF REFORPOST FOR CEMENTATION

1. Position the pin to verify its adaptation. Make a mark on the pin 2 to 3 mm from the antagonist tooth (photo 4) – and cut it with a diamond bit at high rotation under cooling. Do not use scissors or pliers, because the maceration of the fibers changes the mechanical properties of the pin;
  2. Clean the post with alcohol and dry it (photo 5);
  3. Apply a layer of silane (SILANO ANGELUS®) and let it dry for 1 minute (photo 6);
- ① WARNING: It is not necessary to apply the silane on REFORPOST® CARBON FIBER.
4. Apply the ANGELUS® FUSION-DURALINK® CATALYST to the post according to its instructions (photo 7).

#### C. CEMENTATION

1. Clean the canal with the air-water spray and dry it with air and

absorbent paper points (photo 8);

2. Etch the canal with ANGELUS® 37% PHOSPHORIC ACID for 15 seconds (photo 9);
  3. Wash again the canal and dry it with air and absorbent paper points (photo 10). Apply ANGELUS® FUSION-DURALINK® PRIMER in the canal and wait for 1 minute;
  4. Apply ANGELUS® FUSION-DURALINK® CATALYST according to the instructions. The adhesive components must not be mixed before application. They are polymerized chemically without light curing (photo 11);
  5. Spatulate the resinous cement selected\* and fill the conduit;
  6. Put the post in the canal, remove excess of the cement and wait for the polymerization. If the cement is dual, light cure it (photo 13).
- ① Zinc phosphate and glass ionomer cements (conventional and resin-modified) can also be used, but their mechanical properties are inferior to those of the resin cements. Auto-adhesive cements eliminate the needs for the steps of acid conditioning and adhesive application.

#### D. CROWN BUILD-UP

1. Do another etching with ANGELUS® 37% PHOSPHORIC ACID and apply an adhesive on the remaining tooth structures and post. Build the crown core with composite resin (photos 14 and 15), or with ANGELUS® REFORCORE®. REFORCORE® is a prefabricated glass fiber core, for the technique of direct cores (photos 16 and 17).

#### E. USE OF ACCESSORY POSTS IN WEAKENED ROOTS AND WIDE CANALS

1. Fill the remaining spaces of the canal with ANGELUS® REFORPIN® accessory posts. The decrease in thickness of the cement and the better retention provided by REFORPIN® reduce the risk of post/core loosening (photos 18 and 19).

#### WARNINGS | PRECAUTIONS

- Use of direct cores with fiber posts is not indicated for extensive oral rehabilitation in patients with parafunctional habits;
  - Use of the direct core with fiber post should only be indicated to teeth with a minimum height of 2 mm of coronal dentin;
  - Fiber posts should only be cut with water refrigerated high speed diamond burs. Do not use scissors or cutting pliers. Use goggles and a mask to prevent inhalation of particles generated during cutting;
  - Fiber posts must be completely covered by the cementing agent or the composite resin. Exposed fibers absorb water which deteriorates the mechanical properties of the posts.
- ② IMPORTANT: The indication of any post depends on careful evaluation of the quality of endodontic treatment and periapical health.

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**ATTENTION:** This product must be used according to the instructions described in this manual. The manufacturer is not responsible for failure or damage caused by incorrect handling or use.

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For the digital instructions,  
please access this QR code.

GLASS FIBER

For the digital instructions,  
please access this QR code.

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

# ENGLISH - STEEL DRILLS

## **DEFINITION**

STEEL DRILLS are rotary cutting tools used for opening root canals during endodontic treatment and enlarging them to accommodate prosthetic posts.

## **INTENDED APPLICATION**

Preparation of the root canal for endodontic treatment and post cementation.

## **UTILIZATION**

GATES DRILLS: Enlargement of the middle and cervical portions of the root canal.

PEESO DRILLS: Root canal enlargement.

EXACTO DRILLS: For root canal preparation prior to the cementation of EXACTO posts.

## **INTENDED USER**

The steel drills are intended for use by dentists.

## **COMPOSITION**

Stainless steel.

## **PRECAUTIONS**

1. Wear appropriate PPE, such as gloves, eye protection, and surgical masks, to avoid inhaling particles generated by drilling.
2. Before using the steel drill, check for signs of damage or deterioration such as grooves or cracks. Discard any defective drill.
3. Prior to the first and each subsequent use, clean and sterilize the drill. If an instrument or other contaminated product is placed in the sterile container containing the drills, they will become contaminated and must be considered contaminated for future use, storage, transportation, and reuse.
4. Before beginning the procedure, ensure that the equipment (handpiece) is in good working order. Inadequate equipment maintenance may result in user or patient injury.
5. Do not force the steel drill into the equipment. Make sure it is fully seated and locked into the handpiece before using it to avoid damaging the device, patient, or user.
6. Ensure that the drill is moved continuously during use and that no excessive pressure is applied to the drill, as this can cause excessive heat generation and/or damage to the patient and/or user.
7. Only use autoclave sterilization methods because other methods may cause the steel drill to fail prematurely.
8. Do not use corrosive cleaning agents such as chlorine or chloride because the active ingredients are corrosive to stainless steel.
9. Only use the steel drill for the purpose specified by the manufacturer.

## APPLICATION TECHNIQUE

### 1. DRILL INSPECTION

To ensure that the device meets the intended performance and safety standards, it should be inspected before use to identify any damage. When previously unused, the device should be free of cracks, machining drills, machine oils, chatter marks, and device fracture. Particular attention should be paid to the grooves and teeth of the drill. Any device found to have any of the aforementioned flaws should be disposed of in accordance with the instructions provided herein.

The manufacturer recommends that the drills undergo autoclaving a maximum of 5 times. Reuse beyond this has not been validated and is therefore not recommended.

### 2. PRE-TREATMENT

#### Prior to Use

Drills must be inspected for signs of contamination, damage, or deterioration/wear before being used in a dental procedure. Any drill discovered to be in substandard condition should be disposed of in accordance with the guidelines outlined herein.

Steel drills must be sterilized before use according to the instructions below:

#### STERILIZATION PRIOR TO USE

Pack the drills in a surgical grade pouch.

Utilize the following cycle times:

Cycle time	Exposure time (minutes)	Temperature	Drying time (minutes)
PRE-VACUUM (3 PULSES)	3	134°C ± 0	20

#### NOTE:

The instructions from the autoclave manufacturer must be followed and adhered to at all times.

Ensure that the maximum load specified by the sterilizer's manufacturer is not exceeded.

Ensure that the minimum drying time (20 minutes) is not compromised, as failure to do so may result in moisture buildup and corrosion of the drill.

The time between using and cleaning/disinfecting a used steel drill should be less than one hour to reduce the possibility of contaminants drying out and making cleaning difficult. Steel drills should be kept within sterile barriers during the interval between sterilization and use, and only stored and transported with items of the same sterile level. Sterilized drills should never be stored or transported with contaminated items.

### 3. REUSE

**DISINFECTION:** Immerse instruments after use in a special anti-corrosive disinfectant bath and eventually brush them manually.

**RINSING:** Rinse for at least 1 minute.

**CLEANING:** Clean in an ultrasonic bath with detergent.

RINSING: Rinse for at least 1 minute.

DRYING: Dry instruments until there are no visible traces of moisture.

INSPECTION: Carefully inspect each instrument. Sort out damaged or worn-out instruments.

PACKAGING: Place the instruments in sterilization pouches (compliant with en iso 11607-1).

#### Sterilization:

Utilize the following cycle times:

Cycle time	Exposure time (minutes)	Temperature	Drying time (minutes)
Pre-vacuum (3 pulses)	3	134°C ± 0	20

#### NOTE:

The instructions from the autoclave manufacturer must be followed and adhered to at all times.

Ensure that the maximum load specified by the sterilizer's manufacturer is not exceeded.

Ensure that the minimum drying time (20 minutes) is not compromised, as failure to do so may result in moisture buildup and corrosion of the drill.

The time between using and cleaning/disinfecting a used steel drill should be less than one hour to reduce the possibility of contaminants drying out and making cleaning difficult. Steel drills should be kept within sterile barriers during the interval between sterilization and use, and only stored and transported with items of the same sterile level. Sterilized drills should never be stored or transported with contaminated items.

#### Maintenance and inspection:

To ensure proper function and continued safe performance of the drills, inspect them thoroughly after cleaning for any signs of damage and/or deterioration such as corrosion, paying special attention to the grooves and teeth for chips/cracks and the shanks for chatter marks, distortion, and general wear and tear. Any item found in a condition that causes concern must be immediately discarded.

#### Transportation

Steel drills can be transported wet or dry, although if transported wet, there is an increased risk of staining and/or corrosion. In order to prevent damage or deterioration during transportation, suitable protection must be used. Drills must be contained in a clean, dry, and well-maintained drill block/stand or dedicated instrument tray. To reduce the risk of cross contamination, avoid storing clean and soiled drills in the same drill block/stand or instrument tray.

#### Storage

The drills should be stored in the sterilization container (instrument tray, drill stand or pouch) until required. To avoid recontamination of the contents by water, containers or pouches must be dry before opening. Storage should be in dry, clean conditions at room temperature.

Steel drills should be stored in dry, clean containers, such as drill cases and envelopes, that form physical barriers between the drill and the environment when they are kept for extended periods of time, between uses, or during extended periods of non-use.

#### 4. APPLICATION

- a. Place the drill in the hand-piece and tighten it.
- b. Use the procedures specified for each steel drill model, as shown below:
  - GATES DRILLS: Enlargement of the middle and cervical portions of the root canal.
  - PEESO DRILLS: Enlargement of intraradicular conduits.
  - EXACTO DRILLS: For root canal preparation prior to the cementation of EXACTO posts.

#### STORAGE

The product should be stored in a cool, dry, and dirt-free environment.

#### DISPOSAL METHOD

Dispose of steel drills in accordance with national laws and local authority recommendations.

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