RESILAB MASTER

WILCOS

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WILCOS

Product Description

Resilab Master is a highly mimetic light-curing laboratory resin that offers aesthetic and functional characteristics in perfect balance, allowing for the precise reproduction of the original tooth anatomy. It is the first resin in Brazil fully developed according to the Vita Classical shade scale standards. Resilab Master combines advanced technology and modern production processes to achieve highly satisfactory results, meeting the demands of top specialists in Aesthetic Dentistry.

Technical Data:

Average particle size: 0.71 µm (fine-particle hybrid resin) | Vickers hardness: 65.89 MPa | Compressive strength: 168.03 ± 5.56 Kgf | Surface roughness (Ra): 0.141 µm
Radiopacity: Over 50% radiopaque glass | Adhesive strength between Resilab and metals: Gold: 11.26 MPa Cobalt-chromium: 14.77 MPa | Nickel-chromium: 18.70 MPa | Microleakage: Demonstrates excellent marginal integrity with the metals listed above.

Indications: Resilab Master is indicated for the fabrication of anterior and posterior crowns and bridges (applied over a metal framework), inlays, onlays, and laminate veneers.

Composition: Composition: Bis-GMA, Bis-EMA, TEGDMA, UDMA, photoinitiator, fillers, and pigments.

Application Technique: Coating The preparation of coatings must follow the standard stratification method.

Handling: The handling time depends on the ambient temperature and the existing light conditions. Remove the material from the syringe and close it immediately.

Exposure to light can cause premature curing. The approximate handling time for opaque, dentin, enamel, and translucent materials is 30 minutes.

Polymerization: For curing, a light-curing device with a wavelength range of 360 to 500 nm is required, with an optimal maximum of 470 nm.

Layer Thickness: Maximum layer thickness is 1.5 mm per layer, with initial curing (1 minute with Resilab Master) followed by final curing (2 minutes with Resilab Master).

Storage: Store at room temperature (15°C to 30°C), in a dry place, protected from light and heat.

Finishing and Polishing: Finishing and polishing should be performed using small points at a medium speed of 10,000 RPM. The resin surface must first be polished with an appropriate silicone paste. Subsequently, pre-polishing can be performed at 10,000–15,000 RPM using fine pumice powder, a goat hair disc in the handpiece, and a resin polishing paste. Final finishing should be completed with wool discs and a polishing paste for crowns and bridges. Avoid excessive temperatures.

Observations: Whenever possible, work with a suction device. Careful finishing and polishing are essential to achieve optimal results, avoiding residue formation and subsequent color changes.

Precautions, Restrictions, and Warnings:

This product is manufactured for use exclusively as described in the Instructions for Use. Any use outside the described method is the sole responsibility of the professional. This resin is light-sensitive. Excessively intense ambient light may cause premature polymerization. Always keep the product closed when not in use. After being dispensed from its packaging, it should be used promptly. Keep out of reach of children. After use, ensure the material is properly sealed. All dental materials should only be used by licensed dentists and authorized dental laboratories. Restrictions and Warnings There are no contraindications when Resilab Master is used by a dental professional in accordance with the Instructions for Use. However, in patients with a history of allergies to the product's components, its use is not recommended.

Expiration Date: The expiration date listed on the packaging indicates the maximum time for the product to be used. Do not use the product after its expiration date.

Storage: Keep in its original packaging, sealed after use. Store at room temperature (15–30°C) in a dry place, protected from light and heat.

Direct light exposure can alter the properties of Resilab Master. Under proper storage conditions, Resilab Master has a shelf life of 36 months. If stored under refrigeration, the product should be removed and allowed to return to room temperature before use.

RESIBOND

RESIBOND is a primer-type bonding agent used by dental technicians in the metal-resin technique. It cleans the metal surface and creates an adhesion layer between the metal and Resilab Master resin. It is an accessory designed exclusively for use with Resilab Master.

To apply, place a drop of the liquid on a brush or micro-brush and spread it over the metal surface. Allow it to dry and evaporate naturally.

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