Pl

Unimatrix

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En

Unimatrix

The Unimatrix system is indicated for Class II restorations, and composed by preformed sectional matrices and Unimatrix rings. The main advantages of this system offers are the easiness in the reconstruction of contact points and natural tooth anatomy. Unimatrix rings are equipped with silicone tips for a better adaptation and to bring more comfort to the patient.

Use: Isolate the operative field with a rubber dam or other suitable isolation technique. Select, among the 3 matrix sizes, the one that best adapts the cavity to be restored. Install the matrix in the proximal Surface of the tooth to be restored. When it concerns to a MOD cavity, it is advisable to restore one surface at a time. Install the wedge. TDV elastic wedges improve the adaptation of the matrix, in this case, silicone tips are not used and the ring is positioned in the elastic wedge holes. Install the Unimatrix ring. For opening and inserting the ring, TDV Unimatrix pliers or Palfmer pliers can be used. The use or not of the ring and silicone tips, as well as setting up manner depends of the personal preferences and the clinical situation. In concluding the restoration, first remove the ring, then the wedge and finally the matrix. To obtain consistent contact points it is recommended to promote a previous separation of the adjacent tooth, during the initial patient preparation, by inserting a wedge and/or placing the ring.

Sizes: Available in small, medium and large sizes.

Warnings, precautions and restrictions: Professional use. The matrices are for single use only. Their reuse can cause cross contamination. When there is proven allergy to any one of the products components, it must not be used. Excessive expansion with the pliers can lead to deformation of the ring. During the insertion of the ring, maintain the finger pressed on it until the removal of the pliers, ensuring the adjustment and avoiding sliding. Be careful during the ring insertion, due to its flexibility. Make sure that the product is fully seated in the pliers. Avoid leaving the ring open in the pincers for a long period of time. It is recommended to protect the papilla by any method to avoid accidents. To avoid an accidental swallowing or aspiration of this product we always recommend the application of an isolation technique of the operative field. The ring must be used with the triangular silicone tips for an appropriate adjustment of the matrix. Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

RING / SILICONE TIPS - Cleaning: Carry out the cleaning with the aid of water and enzymatic detergent 4 enzymes (Amylase, Protease, Lipase and Carbohydrase) and with pure pH range between 6,5 - 7,5. Prepare the solution according to the manufacturer's instructions. The water should be potable. For a better effectiveness, it is recommended that the tempe-

ture be between 40°C and 55°C. Submerge the products and leave in contact for 5 minutes with the solution. Remove parts and rinse under running water for at least 10 seconds. Dry the pieces with absorbent paper or disposable wipes. **Sterilization:** The ring and the silicone tips must be sterilized separately in autoclave until 126°C - 16 min. The sterilization can cause slight change in the coloring of the tips, but does not affect its functionality. **Inspection:** Check the product before each use. Note if there is damage (breaks/cracks/rust) that prevents proper operation and does not use in such cases. These conditions are indicative of fragility of the material, whose strength can be reduced during clinical use. **Reprocessing Limits:** This is a reusable device. Repeated processing has minimal effect due to its composition. The end of the life of the ring and silicone tips, identified by the professional, is determined due to damage resulting from use.

Composition: Matrices: Stainless steel. Ring: Nickel plated carbon steel. Tubes: Silicone. **Storage:** To be stored in a cool and dry place, protected from sunlight.

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Obelis S.A.