**BLOCO  
PMMA  
MULTI**

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## Inglês

### **Indication:**

*Product indicated for the manufacture of single or multiple crowns, long-lasting screws or cement method, surgical guide, moldable elements, prosthetic bases and others. Product with excellent properties for temporary use lasting.*

### **Instruction for use:**

*The resin body of the block does not require instructions for use, as it is fixed to machining equipment, which brings all the instructions for its operation. There are no restrictions on the equipment that must be used, as it is considered that all equipment on the market is suitable for the intended purposes. This product will no longer be useful if it is not used in CAD/CAM machining equipment.*

*Note: The instructions for using the CAD/CAM machine must be followed according to the manufacturer's instructions.*

### **How to use:**

*The PMMA Block - Zirlux must be handled by a prosthetic professional, who fits the product into a CAD/CAM machining machine, makes adjustments according to the manufacturer's instructions for using such a machine, therefore, the block will be machined with the aim of making teeth for dental prosthesis.*

### **Technical specifications:**

*Zirlux PMMA Blocks are manufactured from the polymerization of methyl methacrylate monomer with poly methyl methacrylate polymer and have the following mechanical properties:*

*Elasticity: 2,220 Mpa*

*Flexural Strength: 100MPa*

*Impact Resistance: 10Kj/m<sup>2</sup>*

*Tensile Strength: 75 Mpa*

*Vickers Curity: 135*

*Water Sorption: 20 micrograms/mm<sup>3</sup>*

*Solubility in water: 0.8microgram/mm<sup>3</sup>*

### **Store:**

*Store away from solvents and intense heat.*

### **Composition:**

*PMMA (Polymethylmethacrylate - Polymer), Ethylene Glycol Dimethacrylate (EGDMA), Fluorescent and Pigments organic.*

*Validity: 10 years after manufacturing.*

*Technical responsible: Rafael Gomes - CRQ: 04263558*

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