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	<b>INSTRUCTION FOR USE</b>	Code: 430020001
	GenDerm	Version / Revision: 5.001 Page: 15 of 42 Publication date: June, 29, 2023

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## INSTRUCTIONS FOR USE

### DESCRIPTION

GenDerm is an implantable medical device, presented in the form of a resorbable bovine collagen membrane to be used by dental and medical professionals in guided bone or tissue regeneration procedures. GenDerm was designed to be used as a biological barrier preventing soft tissue invagination into sites where bone remodeling is taking place. The GenDerm membrane must be hydrated immediately before use (with sterile saline solution) and this makes it flexible and transparent in appearance, and these characteristics favor handling and, consequently, deposition and adaptation over filled areas.

GenDerm membrane is the result of the process of demineralization and lamination of bovine cortical bone. The GenDerm membrane is porous, acellular, biocompatible, non-antigenic, with a high degree of purity, free of contamination by heavy metals. Because it is completely absorbable, the use of GenDerm eliminates the need for a second surgical procedure usually required to remove non-absorbable membranes. In the last production step, the GenDerm membrane goes through the lyophilization process and after packaging it is **STERILIZED** with **GAMMA RADIATION**.

GenDerm is a **SINGLE USE PRODUCT** and **must not be re-sterilized.**

	<b>INSTRUCTION FOR USE</b>	Code: 430020001
	GenDerm	Version / Revision: 5.001 Page: 16 of 42 Publication date: June, 29, 2023

For greater patient comfort, adhesive labels are available inside the final package containing information for product traceability, such as the manufacturing batch number.

## **FUNDAMENTALS OF OPERATION AND ACTION**

One of the goals of periodontal therapy is the correction of anatomical defects or problems caused by the disease process. Several surgical techniques have been suggested in an attempt to regenerate, partially or totally, the tissues lost in advanced periodontal disease. In particular, techniques and tools have evolved to inhibit the growth of adjacent tissues that may harmfully or undesirably replace lost tissues. In recent years, the principle of guided tissue regeneration (GTR) has been developed and improved for the treatment and regeneration of periodontal tissues lost due to inflammatory processes in the periodontium and the regeneration of bone tissue in edentulous areas that would undergo the placement of metallic implants. osseo-integratable. Conceptually, regeneration is a situation in which the repair process, after tissue injury, culminates in the full restoration of the morphological and functional aspects of the tissue (TEN CATE, 1994 apud Taga, 2004).

The principle of RTG is to use a physical barrier to create an adequate environment, protected from the invasion of competing, less differentiated and highly proliferative cells from neighboring tissues, mainly fibroblasts and epithelial cells, and to promote the proliferation of cells originally characteristic of the region for regeneration. from the natural tissues of the area (NYMAN et al., 1980 apud Taga, 2004).

From the 1980s, biocompatible, absorbable and non- absorbable membranes were proposed and began to be commercialized based on the theory that cells originating from the periodontal ligament are the only ones with the potential to differentiate into cementum cells, and the membranes would act as a barrier preventing migration of gingival fibroblasts and epithelial cells over the bare root, allowing selective repopulation by remaining periodontal ligament cells. Currently, the **RTG** technique has also been used in the treatment of angular bone defects around bone-integrated implants.

Two classes of barriers are used in GTR procedures: absorbable and non-absorbable. The indication of which membrane class should be used is directly linked to the type of defect chosen, the cost-benefit ratio and the patient's availability or interest in having a re-entry surgery to remove the non-absorbable membrane. According to what was proposed by Wang et al 1998 apud Taga 2004, such barriers must have

the following characteristics: a) biocompatibility; b) have permeability allowing the diffusion of plasma and nutrients, but not the passage of cells; c) act as a physical support to the adjacent soft tissue, preventing it from collapsing into the space with the blood clot, where cells from the original tissue of the wound can migrate; d) single network protection vascular while the clot undergoes an organization process; e) exclusion of competing cells; f) prevent the formation of a tissue scar; and g) be easy to handle. The advantages of an absorbable material over a non-absorbable one would be:

3- The elimination of the second surgery to remove the barrier, which allows for a reduction in operative time and, consequently, in costs; b) reduction of patient morbidity; c) increased general patient acceptance for RTG and Osteopromotion procedures; and d) reduction of insertion loss risks obtained by regeneration due to re-entry surgery.

4- Biologically absorbable materials have the potential to: a) resist or prevent bacterial colonization; b) increase tissue coverage or protection and reduce barrier exposure; and c) have a greater integration with the tissue host.

GenDerm is a demineralized bovine bone sheet, the organic portion of which is mainly constituted by densely organized type I collagen. The GenDerm production method allows obtaining a biocompatible product with absorption properties by host tissues between 30 and 120 days (depending on the thickness), as has been demonstrated both in animals and in humans.

Collagen membranes are absorbed by phagocytosis and once reabsorbed, they are replaced by endogenous scar tissue and no type of tissue or material foreign to the body is detected. Once hydrated (with sterile saline solution) GenDerm becomes flexible with a transparent appearance, which favors modeling on the graft material, on the bone structure, on the implant or around the implantation region. Once placed, the GenDerm membrane will act as a biological barrier, preventing soft tissue invagination into the implant area and interference with bone regeneration and repair.

Membrane absorption speed can be influenced by local characteristics, such as vascularization of the region, thickness of the mucosa, degree of inflammation resulting from the surgical act itself, changes in the pH of the tissue involved in the repair, local infections and partial exposure to the environment oral.

	<b>INSTRUCTION FOR USE</b>	Code: 430020001
	GenDerm	Version / Revision: 5.001 Page: 18 of 42 Publication date: June, 29, 2023

## CONTENT AND FORMS OF PRESENTATION

GenDerm is presented lyophilized and **STERILE**, individually packaged. The membrane is sold in double blisters: primary blister (rigid plastic heat-sealed in surgical grade paper) packed inside another secondary blister (blister sealed in surgical grade paper) placed in a cardboard box – final packaging. The box is labeled with the data and information needed. It is plasticized, with plastic thermo-shrinkable. The instructions for use of the medical device (instructions for use) are provided inside the final packaging itself – Box.

Traceability – guidelines for healthcare professionals regarding the use and handling of traceability labels.

Self-adhesive labels (five) for product traceability are available inside the final package. The following traceability items are found on the labels: product reference code, manufacturing batch, manufacturer's logo and address, manufacturing date, expiration date, product registration number and other information.

It is obligatory for the dentist or doctor to affix a label (label No. 1) to the patient's clinical record or equivalent document, in order to comply with current health legislation.

To continue the traceability cycle, the manufacturer recommends that the health professional allocate the other available labels (labels n°2 to n°5) as follows:

Etiquette number 2: deliver to patient.


Label number 3: attach to the tax documentation that generates the charge (AIH in the case of a SUS patient, or invoice in the case of a patient treated by the complementary health system).

Label number 4: forward to supplier control.

Label number 5: Forward to the control of the responsible surgeon (main).

It is the entire responsibility of the hospital institution to ensure the identification and traceability of the medical device at the time of use through the adhesive label provided.

	<b>INSTRUCTION FOR USE</b>	Code: 430020001
	GenDerm	Version / Revision: 5.001 Page: 19 of 42 Publication date: June, 29, 2023


GenDerm product as shown.

Outer packaging (left) and membrane inside a double blister (right).

Caption of the graphic symbols contained in the labeling of the medical device.

SYMBOL	MEANING
	Single use product
	Careful, Consult documents escorts
	Fragile, to handle with watch out
	Limit maximum in temperature
	To maintain dry
	Date in manufacturing
	Manufacturer
	Valid until
<b>LOT</b>	Code of Lot
<b>REF</b>	Number of catalog
	Sterilized per Ray Gamma
	Unique sterile barrier system with protective outer packaging
	Do not use it if packaging is damaged.
	Prohibited reesterilization

	<b>INSTRUCTION FOR USE</b>	Code: 430020001
	GenDerm	Version / Revision: 5.001 Page: 20 of 42 Publication date: June, 29, 2023

	Consult Instruction for Use
	Keep Out of the Sun
	Contains biological material of animal origin

Source: ISO 15223 -1 Medical devices - symbols to be used with labels in devices doctors, labeling It is information the be provided - part 1: requirements general

GenDerm membrane is available in the following models, sizes and thicknesses:

**INSTRUCTION FOR USE**

GenDerm

Code: 430020001  
Version / Revision: 5.001  
Page: 21 of 42  
Publication date: June, 29, 2023

GenDerm		
Code	Size	Thickness
980.R.0505	5x5mm	150 - 200 $\mu$
980.R.0506	5 x 6 mm	150 - 200 $\mu$
980.R.0507	5x7mm	150 - 200 $\mu$
980.R.0508	5x8mm	150 - 200 $\mu$
980.R.0509	5x9mm	150 - 200 $\mu$
980.R.0510	5 x 10mm	150 - 200 $\mu$
980.R.1010	10x10 mm	150 - 200 $\mu$
980.R.1015	10 x 15 mm	150 - 200 $\mu$
980.R.1020	10 x 20 mm	150 - 200 $\mu$
980.R.1025	10 x 25 mm	150 - 200 $\mu$
980.R.1030	10 x 30 mm	150 - 200 $\mu$
980.R.1035	10 x 35mm	150 - 200 $\mu$
980.R.1040	10 x 40 mm	150 - 200 $\mu$
980.R.1520	15 x 20 mm	150 - 200 $\mu$
980.R.1525	15 x 25 mm	150 - 200 $\mu$
980.R.1530	15 x 30 mm	150 - 200 $\mu$
980.R.1535	15 x 35mm	150 - 200 $\mu$
980.R.1540	15 x 40 mm	150 - 200 $\mu$
980.S	20 x 20 mm	150 - 200 $\mu$
980.R.2121	21 x 21 mm	150 - 200 $\mu$
980.R.2025	20 x 25 mm	150 - 200 $\mu$
980.R.2030	20 x 30 mm	150 - 200 $\mu$
980.R.2035	20 x 35 mm	150 - 200 $\mu$
980.R.2040	20 x 40 mm	150 - 200 $\mu$
980.R.2530	25 x 30 mm	150 - 200 $\mu$
980.R.2525	25 x 25 mm	150 - 200 $\mu$
980.R.2535	25 x 35 mm	150 - 200 $\mu$
980.R.2540	25 x 40 mm	150 - 200 $\mu$
980.M	30 x 30 mm	150 - 200 $\mu$
980.R.3131	31 x 31 mm	150 - 200 $\mu$
980.R.3035	30 x 35 mm	150 - 200 $\mu$
980.R.3535	35 x 35 mm	150 - 200 $\mu$
980.R.3040	30 x 40 mm	150 - 200 $\mu$
980.R.3540	35 x 40 mm	150 - 200 $\mu$
980.L	40 x 40 mm	150 - 200 $\mu$
980.R.4141	41 x 41mm	150 - 200 $\mu$

GenDerm Flex		
Code	Size	Thickness
983.R.0505	5x5mm	> 200 $\mu$
983.R.0506	5 x 6mm	> 200 $\mu$
983.R.0507	5x7mm	> 200 $\mu$
983.R.0508	5x8mm	> 200 $\mu$
983.R.0509	5x9mm	> 200 $\mu$
983.R.0510	5 x 10mm	> 200 $\mu$
983.R.1010	10x10mm	> 200 $\mu$
983.R.1015	10 x 15 mm	> 200 $\mu$
983.R.1020	10 x 20 mm	> 200 $\mu$
983.R.1025	10 x 25 mm	> 200 $\mu$
983.R.1030	10 x 30 mm	> 200 $\mu$
983.R.1035	10 x 35mm	> 200 $\mu$
983.R.1040	10 x 40 mm	> 200 $\mu$
983.R.1520	15 x 20 mm	> 200 $\mu$
983.R.1525	15 x 25 mm	> 200 $\mu$
983.R.1530	15 x 30 mm	> 200 $\mu$
983.R.1535	15 x 35mm	> 200 $\mu$
983.R.1540	15 x 40 mm	> 200 $\mu$
983.S	20 x 20 mm	> 200 $\mu$
983.R.2121	21 x 21 mm	> 200 $\mu$
983.R.2025	20 x 25 mm	> 200 $\mu$
983.R.2030	20 x 30 mm	> 200 $\mu$
983.R.2035	20 x 35 mm	> 200 $\mu$
983.R.2040	20 x 40 mm	> 200 $\mu$
983.R.2530	25 x 30 mm	> 200 $\mu$
983.R.2525	25 x 25 mm	> 200 $\mu$
983.R.2535	25 x 35 mm	> 200 $\mu$
983.R.2540	25 x 40 mm	> 200 $\mu$
983.M	30 x 30 mm	> 200 $\mu$
983.R.3131	31 x 31 mm	> 200 $\mu$
983.R.3035	30 x 35 mm	> 200 $\mu$
983.R.3535	35 x 35 mm	> 200 $\mu$
983.R.3040	30 x 40 mm	> 200 $\mu$
983.R.3540	35 x 40 mm	> 200 $\mu$
983.L	40 x 40 mm	> 200 $\mu$
983.R.4141	41 x 41mm	> 200 $\mu$

GenDerm Max		
Code	Size	Thickness
984.R.0505	5x5mm	> 300 $\mu$
984.R.0506	5 x 6mm	> 300 $\mu$
984.R.0507	5x7mm	> 300 $\mu$
984.R.0508	5x8mm	> 300 $\mu$
984.R.0509	5x9mm	> 300 $\mu$
984.R.0510	5 x 10 mm	> 300 $\mu$

**Bionext***Hard copy only consultation*

**INSTRUCTION FOR USE**

GenDerm

Code: 430020001  
Version / Revision: 5.001  
Page: 22 of 42  
Publication date: June, 29, 2023

984.R.1010	10x10 mm	> 300 $\mu$
984.R.1015	10 x 15 mm	> 300 $\mu$
984.R.1020	10 x 20 mm	> 300 $\mu$
984.R.1025	10 x 25 mm	> 300 $\mu$
984.R.1030	10 x 30 mm	> 300 $\mu$
984.R.1035	10 x 35mm	> 300 $\mu$
984.R.1040	10 x 40 mm	> 300 $\mu$
984.R.1520	15 x 20 mm	> 300 $\mu$
984.R.1525	15 x 25 mm	> 300 $\mu$
984.R.1530	15 x 30 mm	> 300 $\mu$
984.R.1535	15 x 35mm	> 300 $\mu$
984.R.1540	15 x 40 mm	> 300 $\mu$
984.R.2020	20 x 20 mm	> 300 $\mu$
984.R.2025	20 x 25 mm	> 300 $\mu$
984.R.2030	20 x 30 mm	> 300 $\mu$
984.R.2035	20 x 35 mm	> 300 $\mu$
984.R.2040	20 x 40 mm	> 300 $\mu$
984.R.2530	25 x 30 mm	> 300 $\mu$
984.R.2525	25 x 25 mm	> 300 $\mu$
984.R.2535	25 x 35 mm	> 300 $\mu$
984.R.2540	25 x 40 mm	> 300 $\mu$
984.R.3030	30 x 30 mm	> 300 $\mu$
984.R.3035	30 x 35 mm	> 300 $\mu$
984.R.3535	35 x 35 mm	> 300 $\mu$
984.R.3040	30 x 40 mm	> 300 $\mu$
984.R.3540	35 x 40 mm	> 300 $\mu$
984.R.4040	40 x 40 mm	> 300 $\mu$

**SUPPORT MATERIAL**

The support materials that accompany the medical device are:

- Instructions for use product.
- Traceability and identification labels for protocol and patient.

If the Instruction for Use is not printed with the product, the company provides a two-dimensional barcode (QR Code) on the outside of the package, which can be accessed with reading applications, where it allows access to the Instruction for Use.

**COMPOSITION**

GenDerm is a type I collagen membrane of bovine bone origin. It is an organic material, product of the demineralization of bovine cortical bone until obtaining sheets composed of collagenous protein, as confirmed in X-ray diffraction analysis.

When monitoring the physical-chemical characteristics of GenDerm, the determined components are presented in the following percentages:

Fat = <1%; Protein =>90%; Ash = <1%; pH = 6-8. Heavy metals: As =< 3ppm; Cu = <3ppm; Cd = <5ppm; Cr =<3ppm; Hg = <3ppm; Pb = <30ppm.

**INSTRUCTION FOR USE**

GenDerm

Code: 430020001  
Version / Revision: 5.001  
Page: 23 of 42  
Publication date: June, 29, 2023

**INDICATIONS**

GenDerm is an absorbable, implantable bovine collagen membrane indicated for use in guided tissue regeneration in the covering of other materials in periodontal defects, acting as a biological barrier preventing the invagination of soft tissue to the site where bone repair is taking place. In implant dentistry, it can be used to cover dental implants. This biological barrier can be used in procedures in the areas of orthopedics, neurology and oral and maxillofacial surgery.

**SPECIAL CARE AND CLARIFICATIONS ABOUT THE STORAGE AND TRANSPORTATION OF THE MEDICAL DEVICE**

The packaging must be intact upon receipt.

Open packages must not be used (do not use the medical device if the package is violated);

The product must be transported and stored at room temperature. Avoid humidity and excessive heat;

The medical device storage location must be clean, dry and illuminated in order to maintain ideal storage and transport conditions, as well as its physical and chemical integrity;

The transporter must be informed about the content and deadline of the delivery. The precautions to be taken during transport are shown on the packaging, thus ensuring protection of the medical device from shipment to delivery to the client.

**USE INFORMATION**

Only qualified professionals with knowledge of the surgical techniques and procedures necessary for the proper choice, application and use of the product should use it.

The use of the product with unsuitable surgical techniques and/or procedures and conditions may harm the patient, leading to unsatisfactory results.

The patient must undergo a preoperative evaluation so that the best set of materials is determined and the most appropriate techniques are established for using the product.

postoperative control protocol including the measures taken in the situations encountered.

It is also recommended that the surgeon keep a follow-up report on the evolution of the case together with the results.

In the presence of tumors or infection the use of GenDerm should be avoided.

**INSTRUCTION FOR USE**

GenDerm

Code: 430020001  
Version / Revision: 5.001  
Page: 24 of 42  
Publication date: June, 29, 2023

The adhesive label that accompanies the product must be attached to the patient's protocol to be kept in the clinic or hospital.

Adopt classic prophylactic treatment described in the literature for pre- and immediate and late postoperative periods.

Like all biomaterials, it must be stored in places with low humidity and cool temperatures.

**TECHNICAL AND OPERATIONAL PREPARATION FOR PRODUCT USE**

- The patient must undergo a **preoperative evaluation** so that the best set of materials is determined and the most appropriate techniques are established for using the product.
- The surgeon should adopt **treatment prophylactic** classic described at literature for **pre and immediate postoperative** and late.
- The surgeon should check the **integrity** of the packaging – Do not use it if the packaging is violated.
- The surgeon should check O term in **validity** of product - No use if you are outside of validity.
- The surgical team should be responsible for **preparing** the **patient** and the **area** to perform a **sterile procedure**.
- The surgical team should be responsible for preparing sterile **instruments** and **support materials** (including saline for hydration of the products).
- The team should be responsible for observing **the attire** of the **patient** and the **team** for **sterile procedure**.

**INSTRUCTIONS FOR USE**

- Before opening the package, a container for hydration of the GenDerm membrane must be prepared.
- The **container must be sterile** as well as the fluid used for hydration. It is suggested to use **saline solution sterile**.
- Observe preparation of the area and use of **appropriate sterile clothing** of the doctor, aides and patient for **procedure sterile**.
- Check the integrity and validity of the product before Open
- Using sterile gloves, open the primary blister only in the area prepared for the procedure surgical.
- With a sterile instrument remove the GenDerm membrane from the inside of the Primary package.

**INSTRUCTION FOR USE**

GenDerm

Code: 430020001  
Version / Revision: 5.001  
Page: 25 of 42  
Publication date: June, 29, 2023

- To place for hydrate. To leave at solution until obtaining in aspect transparent and flexible (from 5 to 20 minutes) as shown below.
- With a sterile instrument, position the membrane in place over the implant, over the defect filled.
- Proceed as in custom for closure or coverage from the area implanted.
- A membrane GenDerm he can to be cut in size or Format wanted using scissors sterile.

Detail of the Hydration phase of the GenDerm membrane with saline solution (left). Flexible aspect of the membrane after hydration (right).

**MEDICAL DEVICE DISPOSAL INSTRUCTIONS AND PROCEDURES**

The medical device described in this instruction for use is for **SINGLE USE**, therefore what is not used (leftovers) must be properly disposed of and discarded by the user hospital or dental institution. Leftovers must be autoclaved before final decharacterization by applying a discard cycle or with equivalent lethality (eg 121°C for 30 minutes or 132-134°C 12-10 minutes).

The product whose **packaging** is **damaged** or **inviolate** (compromising the sterility of the product) or with **an expired expiration date, IT SHOULD NOT BE USED.**

The same procedures used to dispose of product leftovers must be followed.

BAUMER recommends that products in the situations described above be mechanically deformed, by making the remains **unusable** by cutting or tearing, after which they must be clearly and visibly identified, indicating their situation, that is, IMPROPER FOR USE.

The methods and procedures for disposing of the used medical device must ensure its **complete de-characterization**, preventing any possibility of its reuse. The mischaracterization of the medical device is the sole responsibility of the medical surgeon or dentist who is using the product, as well as the methods and procedures

	<b>INSTRUCTION FOR USE</b>	Code: 430020001
	GenDerm	Version / Revision: 5.001 Page: 26 of 42 Publication date: June, 29, 2023

for disposal used.

### **CARE AND PRECAUTIONS**

For safe and effective use of the GenDerm Membrane, the following precautions and cautions should be observed:

- Only qualified professionals with knowledge of surgical techniques and procedures necessary for proper selection and application should use the product.
- The medical or dental surgeon should be familiar with the procedures and recommendations for the implantation of GenDerm (see INFORMATION FOR USE).
- Only must to be used instrumentals surgical specific sterilized.
- The surgeon must observe the conditions of preparation and asepsis of the area and instruments, and vestments of the team.
- Long-term steroid therapy, or patients on steroid therapy. anticoagulants.

GenDerm is a **SINGLE USE PRODUCT** and must not be re-sterilized. Opened products must be discarded.

### **WARNINGS**

The surgeon must fully understand not only the medical and surgical aspects of the medical device.

The surgeon must take care to warn or test the patient for sensitivity to products based on bovine derivatives.

Smokers, drug users, alcoholics or malnourished patients should be warned about the increased incidence of bone non-union;

### **POSSIBLE ADVERSE EFFECTS**

Possible complications such as swelling, infections, pain and complications associated with anesthesia may occur as they are common with implant surgeries. As in any type of implant surgery, the patient may experience discomfort in the first postoperative days.

In addition to those mentioned, normal reactions in periodontal surgeries may be found, such as: swelling of the intraoral tissue, thermal sensitivity, gum recession, some loss of crystal bone height, etc. Spontaneous exfoliation can occur if GenDerm has not been properly covered with gingival tissue.

In clinical studies evaluating the GenDerm membrane, no immune reactions or hypersensitivity reactions were found.

**INSTRUCTION FOR USE**

GenDerm

Code: 430020001  
Version / Revision: 5.001  
Page: 27 of 42  
Publication date: June, 29, 2023

**CONTRAINDICATIONS**

The use of GenDerm is not recommended in patients with the following clinical conditions: Uncontrolled metabolism disorders;

Hypersensitivity or allergic reactions to collagen Hypersensitivity or allergic reactions to bovine derivatives Generalized renal dysfunction;

Severe liver disease Bone immaturity;

Osteopenia or severe osteoporosis;

Presence or recent history of systemic or localized infection (contaminated wounds in the area or close to the implant);

Alcoholism or drug abuse;

Pregnancy and lactation (at the discretion of the doctor); fever or leukocytosis;

It is not recommended to use GenDerm in situations where periodontal surgery is not indicated;

The use of GenDerm is not recommended in the presence of tumors or a history of malignancy alone or associated with deproteinized bovine bone matrix (GenOx Inorg) or demineralized (GenOx Org), BMPs (GenPro), composite bone graft (Gen-tech) or other synthetic or natural grafts;

It is not recommended to use GenDerm in any medical or surgical condition that could compromise the success of the surgical procedure.

**STERILITY**

GenDerm medical device is sold **STERILE**. Leftover membranes should not be reused.

**GenDerm is a single-use product and cannot be re-sterilized.**

**GUARANTEE**

The guarantee will apply to the medical device provided it is used in accordance with the provisions of the previous sections of this instruction for use (instruction leaflet): SUPPORT MATERIALS, INDICATIONS, USE INFORMATION, CONTRAINDICATIONS, CAUTIONS AND PRECAUTIONS, RESTRICTIONS, WARNINGS, PACKAGING, CAUTIONS SPECIALS AND CLARIFICATIONS ON THE STORAGE AND TRANSPORT OF THE MEDICAL DEVICE and STERILITY.

**COMPLAINT**

Any customer or user of this medical device who wishes to complain, if they are not satisfied with the services and/or products offered by BAUMER, such as:

**INSTRUCTION FOR USE**

GenDerm

Code: 430020001  
Version / Revision: 5.001  
Page: 28 of 42  
Publication date: June, 29, 2023

Identification, traceability, safety, effectiveness, performance, should contact the authorized BAUMER distributor.

Any defect or suspected defect, or problems in the related medical device, must be immediately reported to BAUMER or an authorized BAUMER distributor.

Any serious adverse effect that affects the health or safety of the patient or user, such as malfunction or damage to the medical device, serious problem or death caused by, or associated with, the use of these products, must be reported to BAUMER.

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