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Biolux OMC is the first two layer acrylic tooth manufactured with sub-nanoparticulated ORMOCER® (ORGanically MODified CERamics) and high-molecular weight Polymethylmethacrylate.

The tooth has high resistance to wear and mastication-impact, chemical attack and color changes caused by alimentation habits.

It is manufacturing within the most modern techniques and under rigid quality control. All production processes comply with ISO 13485:2003 and Good Manufacturing Practices (GMP) norms. They are certified Biocompatible in accordance to ISO 10993-5.

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**SYNTHETIC POLYMER TEETH FOR DENTAL PROSTHESES**

Follow the usual procedures for choice, waxing-up, articulation, polymerization, infection-control and maintenance of the prosthesis, as well as for individual protection.

For an optimized bonding to the denture base resin, grind the tooth base with a grinding instrument. After insulating the flask and prior to pack the acrylic dough, clean the base surface of the tooth using methylmethacrylate monomer.

Indications for the patient: To increase the teeth durability it is highly recommended to avoid contact with alkalis, acids or other substances that can damage acrylic materials during the daily hygiene.

Composition: Polymethylmethacrylate - Ethyleneglycol Dimethacrylate - Organically Modified Ceramics - Fluorescent - Biocompatible Pigments.

Storage: Keep protected from light and heat.

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**P 5** (L5) 30,7 (L6) 8,8 (L7) 33,7 (L8) 8,1

**P 6** (L5) 33,5 (L6) 10,0 (L7) 36,0 (L8) 9,3

Diagram illustrating the measurement points for posterior teeth: L5 (intercusp distance), L6 (tooth width), L7 (intercusp distance), and L8 (tooth width).

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ORMOCER® are substances with the capacity to unite the different materials which are parts of the resin on the base of the same molecule. This avoids the problematic interfaces between the molecules, which come together without forming spaces between them, making the material more resistant.

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