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VITA classical A1 – D4**Instructions for use**

Dear Customer,
Thank you for purchasing the VITA classical A1 – D4 shade guide with
hand-layered ceramic teeth.

The VITA classical A1 – D4 shade guide is used to determine accurate
tooth shade.

Structure

The VITA classical shade family is divided into the following categories:

A1 . A2 . A3 . A3,5 . A4 B1 . B2 . B3 . B4 C1 . C2 . C3 . C4 D2 . D3 . D4
Reddish-brownish reddish-yellowish grayish reddish-grayish

**Procedure**

- Hold the shade guide at a distance of approx. 25-30 cm (arm bent) from
the patient's tooth.
- Decide quickly on the shade by comparing the shade samples of
groups A – D with the natural tooth.
- Disinfect the shade guide after each shade determination.

Tips on determining tooth shade

- Determine tooth shades as much as possible under natural daylight
and daylight lamps and not under typical indoor lighting conditions.
- The surrounding colors should be as neutral as possible. Ask patients to
remove any lipstick or cosmetics which could have an influence on
shade determination. Cover up bright clothing with a gray cape.

– Hold the shade sample tooth as close to the patient's tooth as possible.

– When choosing the shade, only the middle area of the shade sample
tooth should be observed.

– Make your choice quickly and always accept your first decision, since
the eyes begin to tire after approx. 5-7 seconds.

Hygiene tips and instructions for care

The acrylic parts of your VITA classical A1 – D4 shade guide are made of a
high-quality, easy-care material that is temperature resistant (up to
180°C). The entire shade guide can be autoclaved at a max temperature
of 140°C (sterilization according to the gravity method: 30 min at
121°C / 250°F, sterilization according to the pre-vacuum method: 4 min
at 132°C / 270°F).

Disinfection

The VITA classical A1 – D4 shade guide can be surface disinfected. We
recommend customary, mild spray disinfectant or disinfectant wipes.
According to the manufacturers' product information these disinfectants
are effective against HIV and HBV, and are available from most dental
dealers. Please follow the respective manufacturer's instructions for use.
Disinfectants containing phenol, compounds with phenylphenol groups
or methyl ethyl ketone will damage the acrylic parts.

Note

The acrylic parts of the VITA classical A1 – D4 shade guide are neither
monomer, nor UV-resistant, and must be protected from intense
sunlight.

VITA Bleachedguide 3D-MASTER

For planning and monitoring tooth bleaching processes, we recommend
using the VITA Bleachedguide 3D-MASTER. This shade guide was specially
developed for use in the area of tooth bleaching and, as the only shade
guide, offers a gradation of the shade samples in terms of lightness.

VITA Zahnfabrik H. Rauter GmbH & Co.KG
Spitalgasse 3 . D-79713 Bad Säckingen/Germany
Tel. +49 (0) 77 61/5 62-0 . Fax +49 (0) 77 61/5 62-2 99
www.vita-zahnfabrik.com . info @vita-zahnfabrik.com

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