

# BIOSTAR® 4-Aligner

Increase your productivity without compromising on quality.



## A new level in pressure moulding technology.



# A win for your daily laboratory and practice routine!

Productive, precise and flexible - start your own high-end production now, enabling you to thermoform up to 4 aligners at the same time.



## High productivity

Make your workflow more comfortable and efficient. With the BIOSTAR® 4-Aligner it is possible to thermoform up to 4 foils at the same time and thus 160 aligners in one hour.



## Precise results

Each foil is individually heated in its own pressure chamber and thermoformed at 4 bar. The foils fit perfectly to the tooth rim, thus achieving optimum pressure moulding results with the usual high aligner quality.



## Flexible processes

With the BIOSTAR® 4-Aligner, you can choose from various SCHEU-DENTAL foil formats depending on your individual preference: Round foils (Ø 120-125 mm) or rectangular stripes exclusively available from SCHEU-DENTAL for quadruple application (499 x 114 mm).



## Features BIOSTAR® 4-Aligner

BIOSTAR® 4-Aligner	95-230 V
REF	#3030
Power	1450 watt
Working pressure	4 bar
Dimensions (W x H x D)	612 x 337 x 524 mm
Weight	56kg

# An unbeatable combination.

Perfect results with the innovative foil format from SCHEU-DENTAL.



Impressive flexibility: The strips, available exclusively from SCHEU-DENTAL, with a scanable code for easy presetting of the thermoforming parameters, enable you to work efficiently. The BIOSTAR® 4-Aligner sets new standards, especially when using round foils (Ø 120-125 mm): Thanks to the single chamber control, you have the freedom to produce 1, 2, 3 or up to 4 aligners simultaneously as required - customized and economical.

## Best results for different indications:



### Aligner

Material: DURAN®+, CA® Pro+\*, UNIQ pd,  
UNIQ+ pd, IMPRELON® S+ pd, CA® Foil,  
CA® Retention foil pd

\* not available in the USA

**Of course, our service does not end when you purchase one of our devices.**

If you have any questions or experience any problems in your daily working routine, please feel free to contact the SCHEU-DENTAL team. You can call us or arrange an appointment with our field sales staff to discuss specific questions.

