

# Anteriores SUPERIORES

# Anteriores INFERIORES

# POSTERIORES

**LEPTOSSÔMICO**  
Tipo corporal Longilíneo  
Arcada Trapezoidal

	LS	CS	ACS	LCS
2D	42,5	48,0	9,5	7,75
2N	41,6	47,0	10,25	7,5
2P	46,1	52,0	11,25	8,5
A25	44,8	49,5	9,5	8,5
263	42,8	46,0	10,0	7,5
264	45,9	50,0	10,75	8,0
266	49,9	55,0	12,0	9,0

**ATLÉTICO**  
Tipo corporal Intermediário  
Arcada quadrada

	LS	CS	ACS	LCS
A23	39,4	46,0	9,0	7,25
A26	47,3	52,5	9,75	9,0

**PÍCNICO**  
Tipo corporal Arredondado  
Arcada ovalada


	LS	CS	ACS	LCS
3D	44,2	49,0	9,25	8,0
3M	40,0	43,0	9,0	7,0
3N	44,6	48,5	10,0	8,0
3P	48,5	52,5	10,5	8,5

**LI CI ACI LCI**

	LI	CI	ACI	LCI
3M	32,7	33,5	8,25	5,0
2D	32,4	34,0	9,25	5,0
2N	32,7	35,0	9,0	5,0
3D	33,6	36,5	8,75	5,0
3N	33,8	36,5	9,5	5,0
2E	35,4	37,5	9,5	5,5
2P	35,4	38,5	9,75	5,5
46	36,6	39,0	8,75	5,5
3P	37,4	39,5	10,0	6,0
26	37,8	39,5	11,5	5,5

**PSL LMS PIL LMI**

	PSL	LMS	PIL	LMI
30L	30,20	8,3	31,70	9,2
32L	32,05	8,4	33,70	8,6
30M	30,05	8,3	32,15	8,8
32M	32,20	8,3	33,40	8,7
34L	34,15	8,9	35,65	9,4



## Quadro de Articulação

Anteriores Superiores	Anteriores Inferiores	Posteriores Sup. e Inf. 33°
2D	2D	30M - 30L
3D	3D	30M - 32M
2N	2N	30L
2P	2P	34L - 32L
A23	3M	30M
A25	2E	30M - 32M
A26	46	32M
3M	3M	30M
3N	3N	30M - 30L
3P	3P	32M - 32L
263	2N	30M - 30L
264	2E	32M - 32L
266	26	34L

## Tabela de Cores

61	62	65	66	67	69	77	81
----	----	----	----	----	----	----	----

