


W FILE

EMBALAGEM C/ 6 UN

- Tratamento térmico
- Sistema de movimento recíproco
- Simples, segura e eficiente
- Blister estéril, pronto para uso.



	21mm	25mm	31mm
#20	WA1 021 020 BL	WA1 025 020 BL	
#25	WA1 021 025 BL	WA1 025 025 BL	WA1 031 025 BL
#35	WA1 021 035 BL	WA1 025 035 BL	
#45	WA1 021 045 BL	WA1 025 045 BL	
#20-45	WA1 021 A01 BL	WA1 025 A01 BL	

TAPER	#20: 020/7%
	#25: 025/7%
	#35: 035/6%
	#45: 045/5%