



Mixing ratio  
Dosierung  
Rapport de mélange  
Rapporto di miscelazione  
Relación de mezcla  
Proporção de mistura




Mixing time  
Mischungszeit  
Temps de mélange  
Tempo di miscelazione  
Tiempo de mezcla  
Tempo de mistura

**30''**



Working time\*  
Verarbeitungszeit\*  
Temps de travail\*  
Tempo di lavorazione\*  
Tiempo de trabajo\*  
Tempo de processamento\*

**1'30''**



Setting time\*  
Gesamtzeit\*  
Temps de prise\*  
Tempo di indurimento\*  
Tiempo de fraguado\*  
Fim de polimerização\*

**5'30''**



\*23°C (73°F)